





Monopoly Adhesive Bandages

This classic board game is everyone's favorite form of entertainment! Long-lasting adhesive helps bandage stay on, to keep cuts clean. The unique, easy-opening Stat Strip® wrapper allows for aseptic application.

- Perforated, breathable bandage
- Absorbent, non-stick pad for easy removal
- Multi-purpose family protection for cuts and scrapes
- Not made with natural rubber latex
- Sterile



Stat Strip® Opening Technology



- Unique Stat Strip wrapper opening technology
- Patented, easy opening process

Item #	Description	Packaging
1087743	Monopoly, Stat Strip®, 3/4" x 3"	100/bx, 12 bx/cs

HASBRO and its logo, the MONOPOLY name and logo, the distinctive design of the game board, the four corner squares, the MR. MONOPOLY name and character, as well as each of the distinctive elements of the board and playing pieces are trademarks of Hasbro for its property trading game and game equipment. ©1935, 2021 Hasbro. All Rights Reserved.

Licensed By:



Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com