



Better Health. Better Future.

Medical Surgical

Our vast product portfolio provides quality solutions across the healthcare continuum for wherever care is given.



**Enhancing family &
community health
is the way to a
better future**

Traditional Wound Care

Caliber™ Type VII Gauze Sponges	4
Gauze Sponges	4
Gauze Pads	5
Basic Care Gauze Sponges	5
Tray Gauze, Type VII	6
Super Sponges.....	6
X-Ray Detectable Gauze, Type VII	6
Cotton Filled Gauze Sponges.....	6
Clinisorb Non-Woven Sponges	7
Premium Split Sponge.....	7
New Sponges.....	7
Premium Sponges	7
I.V. & T-Drain Sponges.....	8
Top Sponges	8
Basic Care T-Drain Sponges.....	8
Post-Op Sponges.....	8
Eye Pads.....	9
Non-Adherent Pads	9
Pressure Dots.....	9
Non-Adherent Pads with Adhesive.....	9
Emergency Pressure Dressing	10
Trauma Dressing	10
Combine ABD Pads	10
Triangular Bandage	10
Elastic Bandages.....	11
Caliber Non-Woven Island Dressings	11
Premium Elastic Bandages.....	11
Conforming Stretch Gauze.....	12
Basic Care Conforming Stretch Gauze ..	12
Cohesive Bandages	13
Rolled Gauze	13
Antimicrobial Bandage Roll.....	13
Fluff Bandage Roll	14
Basic Care Fluff Bandage Rolls	14
Cotton Tipped Applicators	14
Wood Applicators	14
Tongue Depressors	15
Cotton Balls.....	15
Cotton Roll	15
Cotton Stockinette.....	15
Cotton Undercast Padding	16
Transparent Semi-Permeable Dressings..	16
Caliber Composite Island Dressings.....	16
Medicine Cups	16
Tracheostomy Care Kits.....	17
Suture Removal Kits.....	17

Impregnated Dressings

Petrolatum Gauze.....	17
Packing Strips.....	18
Oil Emulsion	18

Burn Care

Burn Dressings	18
Burn Gel.....	18
Burn Dressings	19

Disposable Blankets & Covers

Ambulance Cot Disposable Linens.....	19
Emergency Blanket.....	20
Survival Wrap	20

Specialty Surgical Sponges

Cherry Sponges.....	20
Kittner Dissectors.....	21
Peanut Sponges.....	21
Prep Sponges.....	21
Flat Stick Sponges	21
Round Stick Sponges	21
Tonsil Sponges.....	21
Vaginal Packing	22
Gauze Packing.....	22

Lap Sponges, OR Towels, & Surgical Drapes

Laparotomy Sponges	22
Caliber™ Surgical Drapes	23
OR Towels.....	23

Diagnostic Products

Lubricating Jelly.....	24
Disposable Vaginal Speculums	24

Tapes & Adhesive Bandages

Transparent Tape.....	24
Paper Tape	24
Cloth Tape	25
Waterproof Tape	25
Caliber™ Single Use Tapes.....	25
Retention Tape	25
Caliber™ Adhesive Bandages	26
Wound Closure Strips	26

Hot and Cold Packs

Instant Cold Pack.....	27
Instant Hot Pack	27

Antiseptics & Cleansers

Alcohol Pads.....	28
BZK Towelettes	28
Cleansing Towelettes.....	28
Nail Polish Remover Pads	28
Sting Relief Pads	29
Lemon Glycerin Swabsticks	29
Povidone Iodine	29

Appendix

HCPCS Codes.....	30-31
Index by Item.....	32-33

Traditional Wound Care

Caliber™ Type VII Gauze Sponges:



- Meets USP Type VII requirements
- 100% woven cotton
- Highly absorbent and offers exceptional clinical performance
- Used for a number of applications including: general wound care, dressing, prepping, wound packing, and debriding
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
C5119	2" x 2", 8-ply	2/pk, 50 pk/bx, 30 bx/cs
C5129	3" x 3", 12-ply	2/pk, 50 pk/bx, 24 bx/cs
C5189	4" x 4", 8-ply	2/pk, 50 pk/bx, 12 bx/cs
C5139	4" x 4", 12-ply	2/pk, 25 pk/bx, 24 bx/cs
C517	8" x 4", 12-ply	2/pk, 25 pk/bx, 20 bx/cs
Non-Sterile		
C82082	2" x 2", 8-ply	200/bg, 25 bg/cs
C82122	2" x 2", 12-ply	200/bg, 40 bg/cs
C83082	3" x 3", 12-ply	200/bg, 20 bg/cs
C84088	4" x 4", 8-ply	200/bg, 20 bg/cs
C84124	4" x 4", 12-ply	200/bg, 10 bg/cs
C84166	4" x 4", 16-ply	200/bg, 10 bg/cs

Gauze Sponges:



- 100% woven cotton
- Ideal for a number of applications including general cleaning, dressing, prepping, packing, and debriding wounds
- Most non-sterile sponges available in 200 or 100 per bag to suit multiple market needs
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
6208	2" x 2", 8-ply	2/pk, 50 pk/bx, 30 bx/cs
6312	3" x 3", 12-ply	2/pk, 50 pk/bx, 24 bx/cs
6408	4" x 4", 8-ply	2/pk, 50 pk/bx, 12 bx/cs
6412	4" x 4", 12-ply	2/pk, 25 pk/bx, 24 bx/cs
6842	8" x 4", 12-ply	2/pk, 25 pk/bx, 20 bx/cs
Non-Sterile		
2085	2" x 2", 8-ply	200/bg, 25 bg/cs
2128	2" x 2", 12-ply	200/bg, 40 bg/cs
2128-100	2" x 2", 12-ply	100/bg, 80 bg/cs
3084	3" x 3", 8-ply	200/bg, 20 bg/cs
3124	3" x 3", 12-ply	200/bg, 20 bg/cs
3124-100	3" x 3", 12-ply	100/bg, 40 bg/cs
4084	4" x 4", 8-ply	200/bg, 20 bg/cs
4122	4" x 4", 12-ply	200/bg, 10 bg/cs
4122-100	4" x 4", 12-ply	100/bg, 20 bg/cs
4162	4" x 4", 16-ply	200/bg, 10 bg/cs
8412	8" x 4", 12-ply	200/bg, 10 bg/cs

Gauze Pads:



- 100% woven cotton
- Individually wrapped for your convenience
- Perforated carton for easy dispensing
- Perfect for when only one pad is needed
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
1212	2" x 2", 12-ply	1/pk, 100 pk/bx, 24 bx/cs
1312	3" x 3", 12-ply	1/pk, 100 pk/bx, 24 bx/cs
1412	4" x 4", 12-ply	1/pk, 100 pk/bx, 12 bx/cs
1212-10	2" x 2", 12-ply	1/pk, 10 pk/bx, 60 bx/cs
1313-10	3" x 3", 12-ply	1/pk, 10 pk/bx, 60 bx/cs
1414-10	4" x 4", 12-ply	1/pk, 10 pk/bx, 60 bx/cs
1212-25	2" x 2", 12-ply	1/pk, 25 pk/bx, 24 bx/cs
1313-25	3" x 3", 12-ply	1/pk, 25 pk/bx, 24 bx/cs
1414-25	4" x 4", 12-ply	1/pk, 25 pk/bx, 24 bx/cs

Basic Care Gauze Sponges:



- 100% woven cotton
- Lighter, economical choice for quality care
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
8500	2" x 2", 8-ply	2/pk, 50 pk/bx, 30 bx/cs
8501	3" x 3", 12-ply	2/pk, 50 pk/bx, 24 bx/cs
8502	4" x 4", 8-ply	2/pk, 50 pk/bx, 12 bx/cs
8503	4" x 4", 12-ply	2/pk, 25 pk/bx, 24 bx/cs
Non-Sterile		
8505	2" x 2", 8-ply	200/bg, 25 bg/cs
8506	2" x 2", 12-ply	200/bg, 40 bg/cs
8507	3" x 3", 8-ply	200/bg, 20 bg/cs
8508	3" x 3", 12-ply	200/bg, 20 bg/cs
8509	4" x 4", 8-ply	200/bg, 20 bg/cs
8510	4" x 4", 12-ply	200/bg, 10 bg/cs

Traditional Wound Care

Tray Gauze, Type VII:



- USP Type VII gauze, 100% woven-cotton
- Packaged in rigid plastic tray with peel back lid
- Sponges are banded in 10's to ensure accurate count
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
408-10	4" x 4", 8-ply	10/ty, 128 ty/cs
412-10	4" x 4", 12-ply	10/ty, 128 ty/cs
416-10	4" x 4", 16-ply	10/ty, 128 ty/cs
421-10	8" x 4", 12-ply	10/ty, 80 ty/cs

Super Sponges:



- Multi-layer, washed, fluff-dried, 100% cotton
- Highly absorbent and durable
- Thick, diagonally folded construction adds loft to enhance overall performance
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
2463	Medium	2/pk, 20 pk/ty, 12 ty/cs
2487	Medium	5/ty, 120 ty/cs
2495	Medium	10/ty, 48 ty/cs
Non-Sterile		
2334	Medium	100/bg, 6 bg/cs

X-Ray Detectable Gauze, Type VII:



- USP Type VII gauze, 100% woven-cotton
- Blue filament runs through each sponge and has a unique wave pattern to ensure easy x-ray detection
- Sterile sponges are packaged in rigid plastic tray with peel back lid
- Both sterile and non-sterile sponges are banded in 10's to ensure accurate count
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
416-2X	4" x 4", 16-ply	10/ty, 128 ty/cs
432-2X	4" x 4", 32-ply	10/ty, 30 ty/cs
482-2X	8" x 4", 12-ply	10/ty, 48 ty/cs
488-2X	8" x 4", 24-ply	10/ty, 48 ty/cs
Non-Sterile		
4412-2X	4" x 4", 12-ply	100/bg, 20 bg/cs
4416-2X	4" x 4", 16-ply	100/bg, 20 bg/cs
4432-2X	4" x 4", 32-ply	50/bg, 20 bg/cs
4812-2X	8" x 4", 12-ply	100/bg, 20 bg/cs
4824-2X	8" x 4", 24-ply	50/bg, 20 bg/cs

Cotton Filled Gauze Sponges:



- 100% woven cotton gauze facing with absorbent cotton fill
- Ideal to help control bleeding following oral surgery or tooth extractions
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
5228	2" x 2", 8-ply	2/pk, 50 pk/bx, 24 bx/cs
2283	2" x 2", 8-ply	200/pk, 25 pk/cs
Non-Sterile		
2280	2" x 2", 8-ply	200/bg, 25 bg/cs

Clinisorb Non-Woven Sponges:



- Rayon/poly blend construction
- Highly absorbent
- Multi-purpose sponges ideal for general wound care and applying ointments
- Virtually non-linting
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
7224	2" x 2", 4-ply	2/pk, 50 pk/bx, 30 bx/cs
7334	3" x 3", 4-ply	2/pk, 40 pk/bx, 30 bx/cs
7444	4" x 4", 4-ply	2/pk, 50 pk/bx, 12 bx/cs
7454SP	4" x 4", 4-ply	5/pk, 20 pk/bx, 8 bx/cs
Non-Sterile		
2102	2" x 2", 4-ply	200/bg, 20 bg/cs
2103	3" x 3", 4-ply	200/bg, 20 bg/cs
2104	4" x 4", 4-ply	200/bg, 10 bg/cs

Premium Split Sponge:



- Unique pre-cut t-slit provides a snug fit around catheters, tubes, I.V.'s and tracheotomies
- Highly absorbent multi-layer (7) sterile dressing protects wounds and reduces dressing changes
- Non-adherent outer layer helps prevent dressing from sticking to wound to reduce the risk of skin trauma
- Thick fluff absorbs fluid and provides cushioning, while cellulose sheets wick away and absorb fluid
- Blue hydrophobic back layer prevents strikethrough, while coloring allows for intuitive dressing placement
- Can be used as a primary or secondary dressing

Item #	Description	Packaging
Sterile		
7420	4" x 6"	1/pk, 10 pk/bx, 12 bx/cs

New Sponges:



- Rayon/poly blend construction
- The effective combination of weight and open aperture allows each sponge to wick more quickly and offer superior absorbency
- Virtually non-linting
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
6122	2" x 2", 4-ply	2/pk, 25 pk/bx, 60 bx/cs
6123	3" x 3", 4-ply	2/pk, 25 pk/bx, 48 bx/cs
6124	4" x 4", 4-ply	2/pk, 25 pk/bx, 24 bx/cs
6125	4" x 4", 4-ply	5/pk, 20 pk/bx, 8 bx/cs
Non-Sterile		
6112	2" x 2", 4-ply	200/bg, 40 bg/cs
6113	3" x 3", 4-ply	200/bg, 20 bg/cs
6114	4" x 4", 4-ply	200/bg, 10 bg/cs

Premium Sponges:



- The most absorbent of our rayon/poly blend sponges
- Heavier fabric offers quick wicking for maximum performance, resulting in a reduction of sponges used
- One of the softest and thickest non-woven sponges on the market
- Ideal for sensitive skin and delicate procedures
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
7224P	2" x 2", 4-ply	2/pk, 50 pk/bx, 30 bx/cs
7334P	3" x 3", 4-ply	2/pk, 40 pk/bx, 30 bx/cs
7444P	4" x 4", 4-ply	2/pk, 50 pk/bx, 12 bx/cs

Traditional Wound Care

I.V. & T-Drain Sponges:



- Pre-cut rayon/poly blend sponges
- Designed to wick moisture away from the site
- Unique pre-cut T-slit provides a snug fit around catheters, tubes, I.V.'s and tracheotomies
- T-Drain sponge also available in 100% woven cotton
- Not made with natural rubber latex

Item #	Description	Packaging
7265	I.V. Sponge, 2" x 2", 6-ply	2/pk, 35 pk/bx, 20 bx/cs
7470	T-Drain Sponge, 4" x 4", 6-ply	2/pk, 25 pk/bx, 12 bx/cs
7468	T-Drain Sponge, 100% Cotton, 4" x 4", 16-ply	1/pk, 25 pk/bx, 12 bx/cs

Basic Care T-Drain Sponges:



- Pre-cut rayon/poly blend sponges
- Lightweight, lower cost alternative
- Designed to wick moisture away from the site
- Unique pre-cut T-slit provides a snug fit around catheters, tubes, I.V.'s and tracheotomies
- Not made with natural rubber latex

Item #	Description	Packaging
8746	4" x 4", 6-ply	2/pk, 25 pk/bx, 12 bx/cs

Post-Op Sponges:



- A multi-layer, woven gauze faced sponge with cellulose fill
- Well suited to be used in place of a single traditional sponge
- Offers fast wicking and increased absorbency
- Not made with natural rubber latex

Item #	Description	Packaging
5434	4" x 3"	2/pk, 25 pk/bx, 24 bx/cs
5444	4" x 4"	2/pk, 25 pk/bx, 24 bx/cs

Top Sponges:



- Gentle non-woven outer facing combined with cellulose liner and rayon fill
- Adds protection to site and offers enhanced absorbency
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
5350	4" x 3"	2/pk, 25 pk/bx, 24 bx/cs
5360	4" x 4"	2/pk, 25 pk/bx, 24 bx/cs
Non-Sterile		
5340	4" x 4"	100/bg, 20 bg/cs

Eye Pads:



- Absorbent cotton filled pad with a soft, gentle non-woven cover that offers protection to this delicate area
- Features an easy place tab on back to ensure clean placement
- Sealed edges provide minimal linting
- Individually wrapped and packaged in a dispenser box for your convenience
- Not made with natural rubber latex

Item #	Description	Packaging
841	1-5/8" x 2-5/8"	1/pk, 50 pk/bx, 12 bx/cs
841B	1-5/8" x 2-5/8"	1/pk, 1000 pk/cs

Non-Adherent Pads:



- Absorbent pad with non-stick outer facing
- Allows dressing to absorb, without adhering to the wound
- Ease of removal helps minimize trauma to the wound site
- Can be cut to fit any size wound without separating
- Not made with natural rubber latex

Item #	Description	Packaging
123	2" x 3"	1/pk, 100 pk/bx, 12 bx/cs
134	3" x 4"	1/pk, 100 pk/bx, 12 bx/cs
138	3" x 8"	1/pk, 75 pk/bx, 8 bx/cs

Pressure Dots:



- Non-woven gauze spheres that provide pressure to insertion site
- An excellent alternative to using a folded 4" x 4" gauze pad
- Not made with natural rubber latex

Item #	Description	Packaging
9550	1" Sphere	1/pouch, 50/bx, 10 bx/cs
9552	1" Sphere	2/pouch, 50/bx, 10 bx/cs

Non-Adherent Pads with Adhesive:



- Highly absorbent cotton pad with a non-adherent film
- Can be used to dress lightly draining or dry sutured wounds, cuts or abrasions
- Adhesive strips keep dressing intact
- Not made with natural rubber latex

Item #	Description	Packaging
124	2" x 3"	1/pk, 100 pk/bx, 24 bx/cs
125	3" x 4"	1/pk, 100 pk/bx, 12 bx/cs

Traditional Wound Care

Emergency Pressure Dressing:



- Convenient all-in-one dressing
- Specifically designed for high absorption and maximum compression
- Second roll can be removed if less pressure is desired
- Not made with natural rubber latex

Item #	Description	Packaging
1050	Emergency Pressure Dressing	1/pk, 50 pk/cs

Trauma Dressing:



- Oversized dressing that offers superior absorbency, protection, and padding to the wound site
- Soft non-woven outer facing wicks moisture into the core which absorbs and disperses fluid quickly to prevent pooling
- Moisture resistant barrier prevents fluid strike-through
- Not made with natural rubber latex

Item #	Description	Packaging
1030TD	10" x 30"	1/pk, 25 pk/cs
3050	12" x 30"	1/pk, 25 pk/cs

Combine ABD Pads:



- Highly absorbent multi-layer dressing
- Soft non-woven outer facing wicks moisture into the cellulose core which quickly absorbs and disperses fluid
- Moisture resistant barrier prevents fluid strike-through
- All edges are sealed to prevent linting and leakage
- Not made with natural rubber latex

Item #	Description	Packaging
5590	5" x 9"	1/pk, 25 pk/ty, 16 ty/cs
5875	8" x 7.5"	1/pk, 20 pk/ty, 12 ty/cs
5810	8" x 10"	1/pk, 20 pk/ty, 16 ty/cs
Non-Sterile		
5940	5" x 9"	25/bg, 16 bg/cs
5942	8" x 7.5"	25/bg, 16 bg/cs
5943	8" x 10"	25/bg, 16 bg/cs
5944	10" x 12"	25/bg, 10 bg/cs
5945	12" x 16"	25/bg, 6 bg/cs
3001	8" x 20 yds Roll	12 rl/cs

Triangular Bandages:



- Makes an excellent all purpose sling
- Natural, woven cotton muslin gauze
- Individually wrapped
- Includes 2 safety pins
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
TB37	37" x 37" x 52"	1/pk, 12 pk/bg, 20 bg/cs
CTB40	40" x 40" x 56"	1/pk, 12 pk/bg, 20 bg/cs

Elastic Bandages:



- Knitted construction offers strength, comfort, and excellent breathability
- Provides firm compression and consistent support
- Individually wrapped
- Includes two metal clips
- Perforated carton for easy dispensing
- LF series not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile Latex		
502	2" x 4.5 yds	10 rl/bx, 5 bx/cs
503	3" x 4.5 yds	10 rl/bx, 5 bx/cs
504	4" x 4.5 yds	10 rl/bx, 5 bx/cs
506	6" x 4.5 yds	10 rl/bx, 5 bx/cs
Non-Sterile		
502LF	2" x 4.5 yds	10 rl/bx, 5 bx/cs
503LF	3" x 4.5 yds	10 rl/bx, 5 bx/cs
504LF	4" x 4.5 yds	10 rl/bx, 5 bx/cs
506LF	6" x 4.5 yds	10 rl/bx, 5 bx/cs

Caliber™ Non-Woven Island Dressings:



- Soft, non-woven backing with adhesive border is gentle on skin
- Highly conformable, breathable, and provides ease of movement
- Absorbent non-adherent pad collects excess exudates and protects wound from trauma

Premium Elastic Bandages:



- Designed with hook and loop closure at both ends, which eliminates the need for clips and provides pressure from the first wrap
- Crocheted construction provides a better hook to secure and minimizes pulls on the fabric for repeat use
- Ideal for situations requiring high compression
- Individually wrapped
- Perforated carton for easy dispensing
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
512LF	2" x 5 yds	10 rl/bx, 5 bx/cs
513LF	3" x 5 yds	10 rl/bx, 5 bx/cs
514LF	4" x 5 yds	10 rl/bx, 5 bx/cs
516LF	6" x 5 yds	10 rl/bx, 5 bx/cs
Sterile		
512LF-S	2" x 5 yds	1/pk, 36 pk/cs
513LF-S	3" x 5 yds	1/pk, 36 pk/cs
514LF-S	4" x 5 yds	1/pk, 36 pk/cs
516LF-S	6" x 5 yds	1/pk, 36 pk/cs

Item #	Description	Packaging
Sterile		
4070	2" x 3"	50/bx, 8 bx/cs
4071	4" x 5"	25/bx, 8 bx/cs
4072	4" x 8"	25/bx, 4 bx/cs
4073	4" x 10"	25/bx, 4 bx/cs
4074	4" x 14"	25/bx, 2 bx/cs
4075	4" x 4"	25/bx, 8 bx/cs
4076	6" x 6"	25/bx, 4 bx/cs

Traditional Wound Care

Conforming Stretch Gauze:



- Rayon/poly knitted stretch gauze bandage
- Provides slight compression to the site while securely holding to any body contour
- Specifically designed not to slip while in use
- Unique weave construction and finished edges prevent unraveling and linting
- Clean wrap is a non-sterile product individually packaged in a polybag
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
701	1" x 4.1 yds	1/pk, 12 pk/bx, 8 bx/cs
702	2" x 4.1 yds	1/pk, 12 pk/bx, 8 bx/cs
703	3" x 4.1 yds	1/pk, 12 pk/bx, 8 bx/cs
704	4" x 4.1 yds	1/pk, 12 pk/bx, 8 bx/cs
706	6" x 4.1 yds	1/pk, 6 pk/bx, 8 bx/cs
Non-Sterile		
601	1" x 4.1 yds	12 rl/bg, 8 bg/cs
602	2" x 4.1 yds	12 rl/bg, 8 bg/cs
603	3" x 4.1 yds	12 rl/bg, 8 bg/cs
604	4" x 4.1 yds	12 rl/bg, 8 bg/cs
606	6" x 4.1 yds	6 rl/bg, 8 bg/cs
602PB	Clean, 2" x 4.1 yds	500/cs
602PB-96	Clean, 2" x 4.1 yds	12 rl/bg, 8 bg/cs
603PB	Clean, 3" x 4.1 yds	500/cs
603PB-96	Clean, 3" x 4.1 yds	12 rl/bg, 8 bg/cs
604PB	Clean, 4" x 4.1 yds	500/cs
604PB-96	Clean, 4" x 4.1 yds	12 rl/bg, 8 bg/cs

Basic Care Conforming Stretch Gauze:



- Rayon/poly knitted stretch gauze bandage
- Provides slight compression to the site while securely holding to any body contour
- Specifically designed not to slip while in use
- Unique weave construction and finished edges prevent unraveling and linting
- Lighter, economical version of the DUKAL Conforming Stretch Gauze
- Single channel construction
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
8513	2"	1/pk, 12 pk/bx, 8 bx/cs
8514	3"	1/pk, 12 pk/bx, 8 bx/cs
8515	4"	1/pk, 12 pk/bx, 8 bx/cs
8516	6"	1/pk, 6 pk/bx, 8 bx/cs
Non-Sterile		
8517	2"	12 rl/bg, 8 bg/cs
8518	3"	12 rl/bg, 8 bg/cs
8519	4"	12 rl/bg, 8 bg/cs
8520	6"	6 rl/bg, 8 bg/cs

Cohesive Bandages:



- Provides consistent support and maintains compression
- Can easily be torn by hand and is self adherent; no clips or fasteners are needed
- Lightweight, porous, and conforms well with movement
- LF series not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile Latex		
8015DB	Dark Blue, 1" x 5 yds	30 rl/bx
8015W	White, 1" x 5 yds	30 rl/bx
8015T	Tan, 1" x 5 yds	30 rl/bx
8015AS	Assorted, 1" x 5 yds	30 rl/bx
8155DB	Dark Blue, 1.5" x 5 yds	48 rl/bx
8155W	White, 1.5" x 5 yds	48 rl/bx
8155T	Tan, 1.5" x 5 yds	48 rl/bx
8155AS	Assorted, 1.5" x 5 yds	48 rl/bx
8025DB	Dark Blue, 2" x 5 yds	36 rl/bx
8025T	Tan, 2" x 5 yds	36 rl/bx
8025AS	Assorted, 2" x 5 yds	36 rl/bx
8035DB	Dark Blue, 3" x 5 yds	24 rl/bx
8035T	Tan, 3" x 5 yds	24 rl/bx
8035AS	Assorted, 3" x 5 yds	24 rl/bx
8045T	Tan, 4" x 5 yds	18 rl/bx
8045AS	Assorted, 4" x 5 yds	18 rl/bx
8065T	Tan, 6" x 5 yds	12 rl/bx
Non-Sterile		
8016TLF	Tan, 1" x 5 yds	30 rl/bx
8016ASLF	Assorted, 1" x 5 yds	30 rl/bx
8156TLF	Tan, 1.5" x 5 yds	48 rl/bx
8156ASLF	Assorted, 1.5" x 5 yds	48 rl/bx
8026TLF	Tan, 2" x 5 yds	36 rl/bx
8026ASLF	Assorted, 2" x 5 yds	36 rl/bx
8036TLF	Tan, 3" x 5 yds	24 rl/bx
8036ASLF	Assorted, 3" x 5 yds	24 rl/bx
8046TLF	Tan, 4" x 5 yds	18 rl/bx
8046ASLF	Assorted, 4" x 5 yds	18 rl/bx
8046BLLF	Blue, 4" x 5 yds	18 rl/bx
8046PLLF	Purple, 4" x 5 yds	18 rl/bx

Rolled Gauze:



- 2-ply 100% woven cotton gauze
- Absorbent, effectively holds dressings in place
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
452	2" x 5 yds	1/pk, 12 pk/bx, 8 bx/cs
453	3" x 5 yds	1/pk, 12 pk/bx, 8 bx/cs
454	4" x 5 yds	1/pk, 12 pk/bx, 8 bx/cs
456	6" x 5 yds	1/pk, 6 pk/bx, 8 bx/cs
Non-Sterile		
402	2" x 5 yds	12 rl/bg, 8 bg/cs
403	3" x 5 yds	12 rl/bg, 8 bg/cs
404	4" x 5 yds	12 rl/bg, 8 bg/cs
406	6" x 5 yds	6 rl/bg, 8 bg/cs

Antimicrobial Bandage Roll



- Contains PHMB and benzalkonium chloride creating a broad spectrum antimicrobial effect within the dressing for gram-positive and gram-negative bacteria and fungi
- May be used as a primary or secondary dressing
- It has been shown to have 4 log bacterial reduction against staphylococcus aureus, escherichia coli, candida albicans, pseudomonas aeruginosa and bacillus subtilis, streptococcus pyogenes, serratia marcescens, aspergillus niger
- As a dressing the fluff roll provides quick wicking action, exceptional absorption, and excellent aeration
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
AB100	4" x 4.1 yd	1/pk, 100 pk/cs

Traditional Wound Care

Fluff Bandage Roll:



- Made of washed, fluff-dried 100% cotton gauze
- Open, crinkle weave provides quick wicking, maximum absorbency, and greater aeration
- Finished edges reduces linting and loose threads
- Adds bulk and loft that will conform as needed, secure dressings in place, and cushion the wound
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
640	3.4" x 3.1 yds, 6-ply	1 rl/bg, 96 rl/cs
645	4.5" x 4.1 yds, 6-ply	1 rl/bg, 100 rl/cs
650	4.5" x 4.1 yds, 8-ply	1 rl/bg, 100 rl/cs
Non-Sterile		
545	4.5" x 4.1 yds, 6-ply	100 rl/cs

Basic Care Fluff Bandage Roll:



- Made of washed, fluff-dried 100% cotton gauze
- Open, crinkle weave provides quick wicking, maximum absorbency, and greater aeration
- Finished edges reduces linting and loose threads
- Adds bulk and loft that will conform as needed, secure dressings in place, and cushion the wound
- Low ply, lightweight alternative
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
8642	2" x 4.5 yds, 3-ply	1 rl/bg, 96 rl/cs
8643	3" x 4.5 yds, 3-ply	1 rl/bg, 96 rl/cs
8645	4.5" x 4.1 yds, 3-ply	1 rl/bg, 100 rl/cs
8644	6" x 4 yds, 3-ply	1 rl/bg, 48 rl/cs

Cotton Tipped Applicators:



- Wood applicators with an absorbent cotton tip
- Sterile 2's available in an easy to peel package
- Non-sterile applicators packaged in bags of 100 with a sterilization indicator strip on the bag
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
9013	3"	2/pk, 100 pk/bx, 10 bx/cs
9016	6"	2/pk, 100 pk/bx, 10 bx/cs
9026	6" - Plastic Shaft	2/pk, 100 pk/bx, 10 bx/cs
Non-Sterile		
9003	3"	100/bg, 10 bg/bx, 10 bx/cs
9006	6"	100/bg, 10 bg/bx, 10 bx/cs

Wood Applicators:



- Smooth 6" applicators
- Excellent for general use
- Packaged in bulk for your convenience
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
9000	6"	1000 pc/bx, 30 bx/cs

Tongue Depressors:



- Clean and smooth hard wood blades
- No odor or taste
- Exceptional for oral use, spreading ointments, mixing medications, or as splints
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
9004	Senior, 6"	1/pk, 100 pk/bx, 10 bx/cs
Non-Sterile		
9001	Junior, 5.5"	500 pc/bx, 10 bx/cs
9002	Senior, 6"	500 pc/bx, 10 bx/cs

Cotton Roll:



- 100% cotton
- Easy to unwind
- Short fibers make it easy to tear
- Packaged in a re-sealable drawstring polybag for your convenience
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
CR1-12	1 lb	1 rl/bg, 12 bg/cs

Cotton Balls:



- 100% cotton
- Soft and uniform in size
- Available in 2 sizes
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
801	Medium	2000/bg, 2 bg/cs
802	Large	1000/bg, 2 bg/cs

Cotton Stockinette:



- 100% natural cotton
- Tubular knitted construction
- Can be used comfortably as a pre-wrap under a cast or to hold bulky dressings in place
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
30112	2" x 25 yds	1 rl/bx, 6 bx/cs
30113	3" x 25 yds	1 rl/bx, 6 bx/cs
30114	4" x 25 yds	1 rl/bx, 6 bx/cs
30116	6" x 25 yds	1 rl/bx, 6 bx/cs

Traditional Wound Care

Cotton Undercast Padding:



- Constructed with 100% soft, cotton fibers
- Provides cushioning between the skin and the cast, while absorbing moisture and allowing the skin to breathe
- Ideal for use with either synthetic or plaster casting material
- Holds padding in place without shifting or bunching
- Easy to tear
- Not made with natural latex

Item #	Description	Packaging
Non-Sterile		
30122	2" x 12'	24 rl/bg, 3 bg/cs
30123	3" x 12'	12 rl/bg, 6 bg/cs
30124	4" x 12'	12 rl/bg, 6 bg/cs
30126	6" x 12'	6 rl/bg, 6 bg/cs

Caliber™ Composite Island Dressings:



- Non-adherent contact layer that covers the wound site
- Soft absorbent pad that collects exudate and protects the wound from trauma
- Non-woven backing with adhesive border that is gentle and holds pad in place
- Waterproof backing to keep wounds dry and keep external contaminants out

Item #	Description	Packaging
4080	4" x 10"	25/bx, 4 bx/cs
4081	4" x 14"	25/bx, 2 bx/cs
4082	4" x 4"	25/bx, 8 bx/cs
4083	6" x 6"	25/bx, 4 bx/cs
4085	4" x 6"	25/bx, 4 bx/cs
4086	6" x 8"	25/bx, 4 bx/cs

Transparent Semi-Permeable Dressings:



- Micro-porous transparent film with high Moisture Vapor Transmission Rate (MVTR)
- Easy two step application
- Protects against fluid contaminants and bacteria
- Rounded edges prevent roll up
- Removable patient label for chart

Item #	Description	Packaging
Sterile		
4045	1-1/2" x 2"	1/pk, 100 pk/bx
4067	2-3/8" x 2-3/4"	1/pk, 100 pk/bx
4112	4" x 4-3/4"	1/pk, 50 pk/bx
4152	6" x 8"	1/pk, 10 pk/bx

Medicine Cups:



- A convenient way to distribute medicine in pill or liquid format
- Made of flexible translucent polypropylene
- Easy to read gradations
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
9007-M	1 oz.	100/slv, 50 slv/cs

Tracheostomy Care Kits:



- Single-use kit complete with all items necessary for care and cleaning of trach site
- Compact and disposable
- Kit offers easy storage, nurse convenience and infection control
- Plastic tray can be used as fluid basin
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
15200	Tracheostomy Care Kit 1 - Pair of Gloves 1 - Trachea Tube Brush 1 - Waterproof Drape 4 - Non-Woven Sponges, 4" x 4" 2 - Pipe Cleaners 2 - Cotton Tipped Applicators, 4" 1 - T-Drain Sponge 1 - Twill Tape, 36"	20 Trays/cs

Suture Removal Kits:



- Single-use kit that includes everything needed for safe and sanitary removal of sutures
- Equipped with littauer scissors, plastic forceps, and gauze pad
- Sterile with an easy to open peel back lid
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
15100	Suture Removal Kit 1 - Littauer Scissors 1 - Plastic Forceps 1 - Gauze Sponge, 3" x 3"	50 Trays/cs

Impregnated Dressings

Petrolatum Gauze:



- Provides soothing non-adherent protection with more convenience and less mess
- Promotes healing by maintaining a moist wound environment
- Clings and conforms to all body contours
- Reduces pain and trauma during dressing changes
- Decreases chances of skin maceration
- Available with Xeroform, a special deodorizing agent containing 3% bismuth tribomphenate that freshens as it protects

Item #	Description	Packaging
Sterile		
204-12	Petrolatum, 1/2" x 72"	12/bx
205-12	Petrolatum, 1" x 36"	12/bx
206	Petrolatum, 3" x 9"	50/bx
207-12	Petrolatum, 3" x 18"	12/bx
208-12	Petrolatum, 3" x 36"	12/bx
209-12	Petrolatum, 6" x 36"	12/bx
210	Petrolatum, 1" x 8"	50/bx
212	Xeroform, 1" x 8"	50/bx
213	Xeroform, 5" x 9"	50/bx
214	Xeroform, 2" x 2"	25/bx
215	Xeroform, 4" x 4"	25/bx
216	Xeroform, 4" x 3 yds	6/bx

Impregnated Dressings

Packing Strips:



- Iodoform provides odor-resistance to the packing strip
- Tamper evident seals
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
260	Plain, 1/4" x 5 yds	12/cs
261	Plain, 1/2" x 5 yds	12/cs
262	Plain, 1" x 5 yds	12/cs
263	Plain, 2" x 5 yds	12/cs
270	Iodoform, 1/4" x 5 yds	12/cs
271	Iodoform, 1/2" x 5 yds	12/cs
272	Iodoform, 1" x 5 yds	12/cs
273	Iodoform, 2" x 5 yds	12/cs

Oil Emulsion:



- Made of a permeable knitted fabric that allows exudates to readily pass through to the second absorbent layer without adhering to granulation tissue
- Enhances wound drainage, supports a wound healing environment, and minimizes maceration
- Flexible and can be cut to size without tearing or unraveling

Item #	Description	Packaging
Sterile		
220	3" x 3"	50/bx
221	3" x 8"	24/bx
222	3" x 8", 3/pk	36/bx
223	3" x 16"	36/bx
224	5" x 9"	12/bx

Burn Care

Burn Dressings:



- Soothes and provides cooling sensation to burns
- Helps prevent wound contamination
- Non-Adherent
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
7007102	4" x 4"	1/pk, 20pk/bx, 12bx/cs
7007103	8" x 8"	1/pk, 50 pk/cs
7007104	16" x 24"	1/pk, 10 pk/cs

Burn Gel:



- Soothes and provides cooling sensation to burns
- Convenient 3.5 gram size foil packets
- Water soluble
- Fragrance free
- Not made with natural rubber latex

Item #	Description	Packaging
7007105-Bulk	Sachet, bulk	1000/cs
7007105-5	Sachet, 5pk	1/pk, 5pk/bg, 200bg/cs
7007105-6	Sachet, 6pk	1/pk, 6pk/bx, 100bx/cs
7007105	Sachet, 25pc box	1/pk, 25pk/bx, 24bx/cs

Burn Dressings:



- 100% woven cotton gauze with thick ply
- Provides superior absorption and a high level of padding and protection to site
- Ideal for debridement of wounds
- Well suited to be used in conjunction with creams or ointments
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
1060	4" x 8", 128-ply	2/pk, 50 pk/cs
9170	18" x 18", 10-ply	1/pk, 100 pk/cs

Disposable Blankets & Covers

Ambulance Cot Disposable Linens:



- All disposable linens are made of durable spunbonded material
- Specifically designed for stretchers and the EMS market
- Sanitary and provides an excellent alternative to the use of traditional linens
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
7100	Pillow Case, Fluid Resistant, 20" x 29", White	25/bg, 4 bg/cs
7102	Flat Sheet, Fluid Resistant, 85" x 41", Blue	5/bg, 20 bg/cs
7103	Flat Sheet, Heavy Duty, Fluid Impervious, 84" x 60", Blue	5/bg, 10 bg/cs
7105	Fitted Sheet, Heavy Duty, Fluid Resistant, 74" x 30" x 22", Blue	5/bg, 20 bg/cs
7106	Fitted Sheet, Extra Heavy Duty, Fluid Resistant, 74" x 30" x 22", Blue	5/bg, 10 bg/cs

Disposable Blankets & Covers

Emergency Blanket:



- Heavy duty, poly coated spunbonded blanket
- Bright yellow for quick identification
- Ideal for protection from the elements in all weather related conditions
- Not made with natural rubber latex

Item #	Description	Packaging
7303	54" x 80"	1/bg, 50 bg/cs

Survival Wrap:



- Heat reflective mylar foil wrap
- Used to trap in heat and maintain body temperature
- Compact packaging allows for convenient storage
- Not made with natural rubber latex

Item #	Description	Packaging
721	52" x 84"	1/pk, 250 pk/cs

Specialty Surgical Sponges

Cherry Sponges:



- 100% woven cotton
- X-ray detectable to ensure visibility
- Used in delicate procedures for dissection and absorption where precision is necessary
- Not made with natural rubber latex

Item #	Description	Packaging
63401-P	1/2", C-5 Holder, 5's	5/pk, 40 pk/cs

Kitner Dissectors and Peanut Sponges:



Kitner Dissectors Peanut Sponges

- 100% woven cotton
- X-ray detectable to ensure visibility
- Ideal for absorption during delicate surgical procedures
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile - Kitner Dissectors		
60404	1/4" x 9/16", No Holder, 5's	5/pk, 40 pk/cs
60401-P	1/4" x 9/16", C-5 Foam Holder, 5's	5/pk, 40 pk/cs
Sterile - Peanut Sponges		
61404	3/8", No Holder, 5's	5/pk, 40 pk/cs
61401-P	3/8", C-5 Foam Holder, 5's	5/pk, 40 pk/cs
Non-Sterile - Peanut Sponges		
61405-5	3/8", C-5 Foam Holder, 5's	5/pk, 50 pk/cs

Prep Sponges:



- 100% woven cotton
- X-ray detectable to ensure visibility
- Highly absorbent
- Can be used to scrub a site with prep solution prior to surgery
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
66406	1-1/4"	100/pk, 5 pk/cs
66407	1-1/2"	100/pk, 5 pk/cs

Flat Stick and Round Stick Sponges:



Flat Stick Sponges

Round Stick Sponges

- 100% woven cotton
- X-ray detectable to ensure visibility
- Highly absorbent
- Ideal as prep sponges prior to surgery
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile - Flat Stick Sponges		
76408	1-1/2"	100/pk, 5 pk/cs
Non-Sterile - Round Stick Sponges		
75405	1", Large, Cotton Fill	100/pk, 5 pk/cs
75406	1-1/2", Super, Cotton Fill	100/pk, 5 pk/cs

Tonsil Sponges:



- 100% woven cotton
- X-ray detectable to ensure visibility
- Used to pack the cavity for absorption
- Double strung for easy retrieval
- Used in tonsillectomies, adenoid surgery, and other procedures
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
74436	Medium, 7/8", on count card	5/pk, 100 pk/cs
74437	Large, 1", on count card	5/pk, 100 pk/cs

Specialty Surgical Sponges

Vaginal Packing:



- 100% woven cotton
- X-ray detectable to ensure visibility
- Best for absorption and packing in cases where there is excessive bleeding
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
1320	4" x 36", 8-ply, Prewashed	1/pk, 100 pk/cs
1339S	2" x 36", 4-ply, Prewashed	1/pk, 100 pk/cs
Non-Sterile		
1336X	4" x 36", 8-ply	200/cs

Gauze Packing:



- 100% woven cotton
- X-ray detectable to ensure visibility
- Absorption and packing for a variety of procedures, wounds, and surgical sites
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
9902	1/2" x 100 yds, 4-ply, Non-Xray	10/cs
9903	1" x 100 yds, 4-ply, Non-Xray	10/cs
9904	2" x 100 yds, 4-ply, Non-Xray	10/cs
9902X	1/2" x 100 yds, 4-ply	10/cs
9903X	1" x 100 yds, 4-ply	10/cs
9904X	2" x 100 yds, 4-ply	10/cs

Lap Sponges, OR Towels, & Surgical Drapes


Laparotomy Sponges:



- Pre-washed, 100% high quality cotton with tight weave
- X-ray detectable
- Banded 5's
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
99-0018	18" x 18"	5/pk, 40 pk/cs
99-0012	12" x 12"	5/pk, 40 pk/cs
99-0036	8" x 36"	5/pk, 40 pk/cs
10-0018	18" x 18"	5/pk, 20 pk/cs
10-0012	12" x 12"	5/pk, 20 pk/cs
10-0004	4" x 18"	5/pk, 20 pk/cs
10-1836-3	18" x 36"	3/pk, 20 pk/cs
Sterile- No Loop, Pre-Treated, CSR Wrap, Soft Pack		
99-9018	18" x 18"	5/pk, 40 pk/cs
Non-Sterile - Pre-Treated		
KP-7018	18" x 18", No Wrap	5/pk, 40 pk/bg, 4 bg/cs
KP-7012	12" x 12"	5/pk, 280 pk/cs

Caliber™ Surgical Drapes:



Item 20-001

- Made of tissue-poly-tissue
- Polyethylene middle layer helps eliminate risk of bacterial migration and moisture strike through
- Available in plain for general use and fenestrated to allow access to operative site while protecting the surrounding area
- Excellent for minor surgical procedures
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
20-001	Fenestrated (3" Round), 18" x 26", 3-ply	1/pk, 50/bx, 6 bx/cs
20-002	Non-Fenestrated, 18" x 26", 3-ply	1/pk, 50/bx, 6 bx/cs

OR Towels:



- Pre-washed, 100% cotton
- Highly absorbent
- All sterile towels are CSR wrapped and have a special corner fold that allows for aseptic presentation into a sterile field
- 17" x 26"
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile - Soft Pouch, CSR Wrapped, Pre-Treated		
CT-01B	Blue	1/pk, 80 pk/cs
CT-02B	Blue	2/pk, 40 pk/cs
CT-04B	Blue	4/pk, 20 pk/cs
CT-04W	White	4/pk, 20 pk/cs
CT-04G	Green	4/pk, 20 pk/cs
CT-06B	Blue	6/pk, 12 pk/cs
CT-08B	Blue	8/pk, 10 pk/cs
CT-12B	Blue	12/pk, 6 pk/cs
Non-Sterile - Bulk		
CT-1730B	Blue	100/cs
Non-Sterile - Bulk, Pre-Treated		
W6020-1	Blue	400/cs

Diagnostic Products

Lubricating Jelly:



- Intended for medical purposes to lubricate body orifices to facilitate entry of diagnostic or therapeutic devices
- Greaseless, water soluble, and bacteriostatic
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
877	2.7g Foil Pack	144 pk/bx, 12 bx/cs
878	5g Foil Pack	150 pk/bx, 10 bx/cs
879	2 oz. tube (flip top)	12 tbs/bx, 6 bx/cs
880	4 oz. tube (flip top)	12 tbs/bx, 6 bx/cs
881	1 oz. tube (flip top)	12 tbs/bx, 6 bx/cs

Disposable Vaginal Speculums:




- Made of heavy grade plastic with smooth edges and offers a wide range of adjustability
- Disposable; helps reduce the risk of cross contamination
- Graves style
- All speculums are individually wrapped for hygienic purposes
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
6650	Small	10/bg, 10 bg/cs
6660	Medium	10/bg, 10 bg/cs
6670	Large	10/bg, 10 bg/cs

Tapes & Adhesive Bandages

Transparent Tape:



- High quality porous clear tape that is perforated to be easily torn bi-directionally
- Excellent for I.V. sites, dressings, and tubing retention
- Conforms and stretches to accommodate movement
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
T50	1/2" x 10 yds	24 rl/bx, 12 bx/cs
T110	1" x 10 yds	12 rl/bx, 12 bx/cs
T210	2" x 10 yds	6 rl/bx, 12 bx/cs
T310	3" x 10 yds	4 rl/bx, 12 bx/cs

Paper Tape:



- Light-weight, microporous paper tape that will hold even when wet
- Gentle on skin and won't harm delicate wounds
- Recommended for delicate and elderly skin
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
P50	1/2" x 10 yds	24 rl/bx, 12 bx/cs
P110	1" x 10 yds	12 rl/bx, 12 bx/cs
P210	2" x 10 yds	6 rl/bx, 12 bx/cs
P310	3" x 10 yds	4 rl/bx, 12 bx/cs

Cloth Tape:



- Durable cloth tape designed for strength
- Gentle adhesive will leave minimal residue on the skin
- Designed for bilateral tearing without use of scissors
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
C50	1/2" x 10 yds	24 rl/bx, 12 bx/cs
C110	1" x 10 yds	12 rl/bx, 12 bx/cs
C210	2" x 10 yds	6 rl/bx, 12 bx/cs
C310	3" x 10 yds	4 rl/bx, 12 bx/cs

Waterproof Tape:



- Provides high tensile strength with a secure hold
- Offered in plastic spools and shells
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
AT2316	1/2" x 5 yds	24 rl/bx, 18 bx/cs
AT2323	1/2" x 2.5 yds	36 rl/bx, 16 bx/cs
AT110	1" x 10 yds, bulk	144/cs
AT150	1" x 5 yds, bulk	216/cs

Caliber™ Single Use Tapes:



- Our high quality tapes in patient specific, single use rolls
- The perfect amount of tape for securing dressings, I.V. sites, or tubing
- Contributes to the safe and sanitary treatment of all patients
- Available in paper, cloth, or transparent
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
P15	Paper, 1" x 1.5 yds	100 rl/bx, 5 bx/cs
T15	Transparent, 1" x 1.5 yds	100 rl/bx, 5 bx/cs
C15	Cloth, 1" x 1.5 yds	100 rl/bx, 5 bx/cs

Retention Tape:



- Soft & flexible making it ideal for difficult to tape areas
- Printed s-curve release liner
- Pre-cut to eliminate waste
- Provides light compression without constriction
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
8220	2" x 11 yds	1 rl/bx, 10 bx/cs
8221	4" x 11 yds	1 rl/bx, 6 bx/cs
8222	6" x 11 yds	1 rl/bx, 4 bx/cs

Tapes & Adhesive Bandages

Caliber™ Adhesive Bandages:




- Durable plastic, sheer, or flexible fabric to meet all market needs
- Absorbent, non-adherent pad protects and cushions the wound without sticking
- Sterile and individually wrapped in easy to peel pouches
- Not made with natural rubber latex



Item #	Description	Packaging
Sterile - Plastic		
7616	3/4" x 3"	100/bx, 24 bx/cs
7617	1" x 3"	100/bx, 24 bx/cs
7611	2" x 4", X-Large	100/bx, 24 bx/cs
7614	Spot, 7/8" Round	100/bx, 24 bx/cs
Sterile - Clear/Sheer		
7604	3/4" x 3"	100/bx, 24 bx/cs
7608	1" x 3"	100/bx, 24 bx/cs
7612	2" x 4", X-Large	100/bx, 24 bx/cs
7615	Spot, 7/8" Round	100/bx, 24 bx/cs
Sterile - Fabric		
7602	Knuckle	100/bx, 24 bx/cs
7605	3/4" x 3"	100/bx, 24 bx/cs
7609	1" x 3"	100/bx, 24 bx/cs
7613	2" x 4", X-Large	100/bx, 24 bx/cs
7618	Fingertip	100/bx, 24 bx/cs
Sterile - Butterfly Wound Closure Strips		
7600	3/8" x 1-13/16"	100/bx, 24 bx/cs

Wound Closure Strips:



- Adhesive backed, thin, non-woven, reinforced fibers that are arranged in a bilateral pattern for high tensile strength
- Provides support following early suture/staple removal
- Keeps wound closed without significant tension, edema, stretch or movement
- Translucent and porous
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
5150	1/8" x 3"	5/pk, 50 pk/bx, 4 bx/cs
5151	1/4" x 3"	3/pk, 50 pk/bx, 4 bx/cs
5152	1/4" x 1-1/2"	6/pk, 50 pk/bx, 4 bx/cs
5156	1/4" x 4"	10/pk, 50 pk/bx, 4 bx/cs
5157	1/2" x 4"	6/pk, 50 pk/bx, 4 bx/cs
5158	1" x 4"	4/pk, 50 pk/bx, 4 bx/cs
Non-Sterile		
C-02	1/4" x 3"	500/bx
C-03	1/4" x 1-1/2"	500/bx
C-04	1/4" x 4"	500/bx
C-05	1/2" x 4"	500/bx
C-06	1" x 4"	500/bx

Instant Cold Pack:



- Flexible material conforms to body
- Helps reduce swelling
- Ideal for minor sprains, muscle spasms, headaches, and burns

Item #	Description	Packaging
9850	6" x 8"	24/cs
9851	4" x 5", Junior	50/cs

Instant Hot Pack:



- Effective pain relief
- Consistent heat for up to 30 minutes
- Activates immediately
- Flexible pack conforms to body

Item #	Description	Packaging
9840	5" x 9"	24/cs

Antiseptics & Cleansers

Alcohol Pads:



- Saturated with 70% isopropyl alcohol
- Ideal for use as an antiseptic cleanser
- 2-ply, non-woven pad
- Not made with natural rubber latex

Item #	Description	Packaging
Sterile		
853	Medium	1/pk, 200/bx, 20 bx/cs
857-4000	Medium, Bulk	1/pk, 4000/cs
861	Large	1/pk, 100/bx, 10 bx/cs
Non-Sterile		
852	Medium	1/pk, 200/bx, 20 bx/cs
852-1000	Medium, Bulk	1/pk, 1000/bx, 4 bx/cs

BZK Towelettes:



- Helps to reduce the risk of infection in minor cuts, scrapes and burns
- 5" x 8"
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
855	5" x 8"	1/pk, 100/bx, 20 bx/cs
855-1000	5" x 8", Bulk	1/pk, 1000/cs
855-1200	5" x 8", Bulk	3/strip, 1/pk, 1200/cs

Cleansing Towelettes:



- For cleansing of face, hands, and body without soap and water
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
858	5" x 8"	1/pk, 100/bx, 20 bx/cs
858-1000	5" x 8", Bulk	3/strip, 1/pk, 1500/cs

Nail Polish Remover Pads:



- Effectively removes nail polish
- 2-ply non-woven pad
- Acetone Free
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
862	Medium	1/pk, 100/bx, 10 bx/cs

Sting Relief Pads:



- Temporary relief of pain and itching associated with minor burns, scrapes, and insect bites
- 2-ply non-woven pad
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
856	Medium	1/pk, 200/bx, 20 bx/cs
856-1000	Medium, Bulk	3/strip, 1/pk, 1500/bx, 2 bx/cs

Lemon Glycerin Swabsticks:



- Pre-moistened with refreshing lemon flavored glycerin
- For post operative oral care, relief of dry mouth, and swabbing when special mouth care is prescribed
- Triples
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
864	Swabsticks, 4"	3/pk, 25 pk/bx, 10 bx/cs

Povidone Iodine:



- Povidone-Iodine USP 10% solution
- Used for patient preoperative skin preparation
- Helps to reduce bacteria that potentially can cause skin infection
- Sterile, single-use, easy to use foil pouch

Item #	Description	Packaging
Sterile		
882-4000	PVP-I Prep Pads, Medium	4000/cs
883-500	PVP-I Prep Swabstick 3's, 4"	500/cs
884-1000	PVP-I Prep Swabstick 1's, 4"	1000/cs
885-500	PVP-I 3/4 oz Solution Pouch	500/cs
886-250	PVP-I 1 oz Solution Pouch	250/cs

Item #	HCPCS	Item #	HCPCS	Item #	HCPCS	Item #	HCPCS
123	A6251	645	A6446	4074	A6254	6125	A6402
124	A6254	650	A6446	4075	A6254	6208	A6402
125	A6254	701	A6445	4076	A6254	6312	A6402
134	A6251	702	A6445	4080	A6203	6408	A6402
138	A6252	703	A6446	4081	A6204	6412	A6402
206	A6223	704	A6446	4082	A6203	6842	A6403
210	A6222	706	A6447	4083	A6203	7224	A6402
212	A6222	841	A6410	4084	A6216	7265	A6402
213	A6223	852	A4245	4085	A6203	7334	A6402
220	A6222	853	A4245	4086	A6204	7420	A6252
221	A6223	861	A4245	4112	A6258	7444	A6402
222	A6223	877	A4332	4122	A6216	7468	A6402
223	A6223	878	A4332	4152	A6258	7470	A6402
260	A6407	879	A4402	4162	A6216	7600	A4450
261	A6407	880	A4402	5150	A4450	7602	A6413
262	A6407	1050	A6446	5151	A4450	7604	A6413
263	A6407	1060	A6403	5152	A4450	7605	A6413
270	A6266	1212	A6402	5156	A4450	7608	A6413
271	A6266	1312	A6402	5157	A4450	7609	A6413
272	A6266	1412	A6402	5158	A4450	7611	A6413
273	A6266	2085	A6216	5228	A6402	7612	A6413
402	A6442	2102	A6216	5340	A6216	7613	A6413
403	A6442	2103	A6216	5350	A6402	7614	A6413
404	A6443	2104	A6216	5360	A6402	7615	A6413
406	A6444	2128	A6216	5434	A6402	7616	A6413
452	A6445	2280	A6216	5444	A6402	7617	A6413
453	A6445	2283	A6402	5590	A6252	7618	A6413
454	A6446	2334	A6216	5810	A6253	8220	A4450
456	A6447	2463	A6402	5875	A6252	8221	A4450
502	A6448	2487	A6402	5940	A6217	8222	A4450
503	A6448	2495	A6402	5942	A6218	8412	A6217
504	A6449	3050	A6404	5943	A6218	8500	A6402
506	A6450	3084	A6216	5944	A6218	8501	A6402
545	A6443	3124	A6216	5945	A6218	8502	A6402
601	A6442	4045	A6257	6112	A6216	8503	A6402
602	A6442	4067	A6257	6113	A6216	8505	A6216
603	A6442	4070	A6219	6114	A6216	8506	A6216
604	A6443	4071	A6254	6122	A6402	8507	A6216
606	A6444	4072	A6254	6123	A6402	8508	A6216
640	A6446	4073	A6254	6124	A6402	8509	A6216

Item #	HCPCS	Item #	HCPCS	Item #	HCPCS	Item #	HCPCS
8510	A6216	207-12	A6224	8015T	A6453	C83082	A6216
8513	A6445	208-12	A6224	8015W	A6453	C84088	A6216
8514	A6446	209-12	A6224	8016ASLF	A6453	C84124	A6216
8515	A6446	2128-100	A6216	8016TLF	A6453	C84166	A6216
8516	A6447	3124-100	A6216	8025AS	A6453	KP-7012	A6218
8517	A6442	408-10	A6402	8025DB	A6453	KP-7018	A6218
8518	A6443	412-10	A6402	8025T	A6453	P110	A4450
8519	A6443	4122-100	A6216	8026ASLF	A6453	P15	A4450
8520	A6444	416-10	A6402	8026TLF	A6453	P210	A4450
8642	A6445	416-2X	A6402	8035AS	A6453	P310	A4450
8643	A6446	421-10	A6403	8035DB	A6453	P50	A4450
8644	A6447	432-2X	A6402	8035T	A6453	T110	A4450
8645	A6446	4412-2X	A6216	8036ASLF	A6453	T15	A4450
8746	A6402	4416-2X	A6216	8036TLF	A6453	T210	A4450
9170	A6404	4432-2X	A6216	8045AS	A6454	T310	A4450
11724	A4554	4812-2X	A6217	8045T	A6454	T50	A4450
12424	A4554	482-2X	A6403	8046ASLF	A6454		
12436	A4554	4824-2X	A6217	8046TLF	A6454		
15200	A4629	488-2X	A6403	8155AS	A6453		
30112	A6457	502LF	A6448	8155DB	A6453		
30113	A6457	503LF	A6449	8155T	A6453		
30114	A6457	503LF-B	A6449	8155W	A6453		
30116	A6457	504LF	A6449	8156ASLF	A6453		
30122	A9270	506LF	A6450	8156TLF	A6453		
30123	A9270	512LF	A6448	841B	A6410		
30124	A9270	513LF	A6449	852-1000	A4245		
30126	A9270	514LF	A6449	AT2316	A4452		
81724	A4554	516LF	A6450	AT2323	A4452		
1030TD	A6404	602PB	A6442	C110	A4450		
1212-10	A6402	602PB-96	A6442	C15	A4450		
1212-25	A6402	603PB	A6442	C210	A4450		
1212-36B	A6402	603PB-96	A6442	C310	A4450		
1312-36B	A6402	604PB	A6443	C50	A4450		
1313-10	A6402	604PB-96	A6443	C5119	A6402		
1313-25	A6402	7224P	A6402	C5129	A6402		
1412-18B	A6402	7334P	A6402	C5139	A6402		
1414-10	A6402	7444P	A6402	C517	A6403		
1414-25	A6402	7454SP	A6402	C5189	A6402		
204-12	A6223	8015AS	A6453	C82082	A6216		
205-12	A6223	8015DB	A6453	C82122	A6216		

Index

Item #	Page	Item #	Page	Item #	Page	Item #	Page	Item #	Page
123	9	453	13	1412	5	5152	26	7420	7
124	9	454	13	1540	30	5156	26	7444	7
125	9	456	13	1541	30	5157	26	7468	8
134	9	502	11	1542	30	5158	26	7470	8
138	9	503	11	1550	30	5228	6	7600	26
206	17	504	11	1560	30	5340	8	7602	26
210	17	506	11	2085	4	5350	8	7604	26
212	17	545	14	2102	7	5360	8	7605	26
213	17	601	12	2103	7	5434	8	7608	26
214	17	602	12	2104	7	5444	8	7609	26
215	17	603	12	2128	4	5590	10	7611	26
216	17	604	12	2280	6	5810	10	7612	26
220	18	606	12	2283	6	5875	10	7613	26
221	18	640	14	2334	6	5940	10	7614	26
222	18	645	14	2463	6	5942	10	7615	26
223	18	650	14	2487	6	5943	10	7616	26
224	18	701	12	2495	6	5944	10	7617	26
260	18	702	12	3001	10	5945	10	7618	26
261	18	703	12	3050	10	6112	7	8220	25
262	18	704	12	3084	4	6113	7	8221	25
263	18	706	12	3124	4	6114	7	8222	25
270	18	721	20	4045	16	6122	7	8412	4
271	18	801	15	4067	16	6123	7	8500	5
272	18	802	15	4070	11	6124	7	8501	5
273	18	841	9	4071	11	6125	7	8502	5
301	29	852	32	4072	11	6208	4	8503	5
303	29	853	32	4073	11	6312	4	8505	5
304	28	855	32	4074	11	6408	4	8506	5
305	28	856	34	4075	11	6412	4	8507	5
311	31	858	33	4076	11	6650	24	8508	5
340	29	861	33	4080	16	6660	24	8509	5
341	29	862	33	4081	16	6670	24	8510	5
342	29	864	34	4082	16	6842	4	8513	12
343	29	877	24	4083	16	7100	19	8514	12
344	29	878	24	4085	16	7102	19	8515	12
350	31	879	24	4086	16	7103	19	8516	12
352	31	880	24	4084	4	7105	19	8517	12
360	28	881	24	4112	16	7106	19	8518	12
402	12	1050	10	4122	4	7107	19	8519	12
403	12	1060	19	4152	16	7224	7	8520	12
404	12	1212	5	4162	4	7265	8	8642	14
406	12	1312	5	5150	26	7303	20	8643	14
452	13	1320	22	5151	26	7334	7	8644	14

Item #	Page	Item #	Page	Item #	Page	Item #	Page	Item #	Page
8645	14	10-0004	22	514LF	11	8046BLLF	13	C310	25
8746	8	10-0012	22	516LF	11	8046PLLF	13	C50	25
9000	14	10-0018	22	512LF-S	11	8046TLF	13	C5119	4
9001	15	10-1836-3	22	513LF-S	11	8065T	13	C5129	4
9002	15	1030TD	10	514LF-S	11	8155AS	13	C5139	4
9003	14	1212-10	5	516LF-S	11	8155DB	13	C517	4
9004	15	1212-25	5	602PB	12	8155T	13	C5189	4
9006	14	1313-10	5	602PB-96	12	8155W	13	C82082	4
9013	14	1313-25	5	603PB	12	8156ASLF	13	C82122	4
9016	14	1336X	22	603PB-96	12	8156TLF	13	C83082	4
9026	14	1339S	22	60401-P	21	841B	9	C84088	4
9170	19	1414-10	5	604PB	12	852-1000	33	C84124	4
9550	9	1414-25	5	604PB-96	12	855-1000	33	C84166	4
9552	9	204-12	17	61401-P	21	855-1200	33	CR1-12	15
9850	32	205-12	17	61405-5	21	856-1000	34	CT-01B	23
9851	32	207-12	17	63401-P	20	857-4000	28	CT-02B	23
9840	32	208-12	17	7224P	8	858-1000	33	CT-04B	23
9902	22	209-12	17	7334P	8	882-4000	29	CT-04G	23
9903	22	20-001	23	7444P	8	883-500	29	CT-04W	23
9904	22	20-002	23	7454SP	7	884-1000	29	CT-06B	23
15100	17	2128-100	4	7007105 - Bulk	18	885-500	29	CT-08B	23
15200	17	3124-100	4	7007105-5	18	886-250	29	CT-12B	23
30112	15	408-10	6	7007105-6	18	9007-M	16	CT-1730B	23
30113	15	412-10	6	8015AS	13	99-0012	22	CTB40	10
30114	15	4122-100	4	8015DB	13	99-0018	22	KP-7012	22
30116	15	416-10	6	8015T	13	99-0036	22	KP-7018	22
30122	16	416-2X	6	8015W	13	9902X	22	P110	24
30123	16	421-10	6	8016ASLF	13	9903X	22	P15	24
30124	16	432-2X	6	8016TLF	13	9904X	22	P210	24
30126	16	4412-2X	6	8025AS	13	99-9018	23	P310	24
60404	20	4416-2X	6	8025DB	13	AB100	13	P50	24
61404	20	4432-2X	6	8025T	13	AT110	25	T110	24
66406	20	4812-2X	6	8026ASLF	13	AT150	25	T15	24
66407	20	482-2X	6	8026TLF	13	AT2316	25	T210	24
74436	21	4824-2X	6	8035AS	13	AT2323	25	T310	24
74437	21	486-2X	6	8035DB	13	C-02	26	T50	24
75405	21	488-2X	6	8035T	13	C-03	26	TB37	10
75406	21	502LF	11	8036ASLF	13	C-04	26	W6020-1	24
76408	21	503LF	11	8036TLF	13	C-05	26		
7007102	18	504LF	11	8045AS	13	C-06	26		
7007103	18	506LF	11	8045T	13	C110	25		
7007104	18	512LF	11	8046ASLF	13	C15	25		
7007105	18	513LF	11			C210	25		



Dukal, LLC

2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800

Toll Free: 1 800 243 0741

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Find our full line of products at [dukal.com](https://www.dukal.com)