



Better Health. Better Future.

Pandemic Preparedness

Products that provide protection for healthcare professionals working in infectious conditions.



Enhancing family & community health is the way to a better future

Gowns

- Isolation Gowns4
- CPE Gowns6
- AAMI Level 1 Gowns.....8
- AAMI Level 2 Gowns.....10

Disposable Scrubs and Coveralls

- Coversoff™ Disposable Apparel.....12
- Disposable Scrub Wear13
- Disposable Coveralls14

Caps and Shoe Covers

- Bouffant Caps.....15
- Shoe Covers.....16

Lab Coats

- Lab Coats.....17

Patient Exam Wear

- Patient Exam Gowns.....18
- Patient Exam Shorts.....18

Face Masks

- General Use Face Mask19
- Surgical Masks20
- Face Shield22

Diagnostics and Instruments

- Digital Thermometer23
- Disposable Stethoscope24
- Bowles Stethoscope.....25
- Single Head Stethoscope26
- Dual Head Stethoscope.....27
- Liquid Dispenser.....28
- Non-Latex Tourniquet29
- Quick Release Tourniquet30
- Finger Cots31
- Nitrile Finger Cots.....32
- IV Stands.....33

Post Mortem

- Body Bag.....34



Isolation Gowns

Made of high quality spunbonded polypropylene

Style Guide



Wrist Style



Neck Style



- Each gown features elastic cuffs with waist and neck tie closures
- Available in multiple colors to allow easy detection between departments or functions



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
301	Isolation Gown, Yellow	10/bg, 5 bg/cs
301BL	Isolation Gown, Blue	10/bg, 5 bg/cs

Poly Coated Isolation Gowns

Made of high quality spunbond polypropylene with a polyethylene coating

Style Guide



Wrist Style



Neck Style



- Each gown features elastic cuffs with waist and neck tie closures
- Available in multiple colors to allow easy detection between departments or functions



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
303	Isolation Gown, coated, Yellow	10/bg, 5 bg/cs
303BL	Isolation Gown, coated, Blue	10/bg, 5 bg/cs

CPE Isolation Gowns

Polyethylene film provides enhanced fluid protection

- Over-the-head protective cover gown with thumb loops
- Tear away design enables gowns to be removed quickly
- Thumb loop holes provide easy glove doning
- Open back style provides comfortable fit
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
307R	CPE Gown, Regular	15/bg, 5 bg/cs
307XL	CPE Gown, XL	15/bg, 5 bg/cs

Style Guide



Wrist Style



Neck Style



AAMI LEVEL 1

Isolation Gowns

Made of multi-ply SMS material

- Meets AAMI Level 1 requirements for fluid protection
- Provides a slight barrier against small amounts of fluid penetration
- Elastic wrist with thumb loops to provide a more secure fit under gloves
- Tape tab neck with waist ties
- Used in minimal risk situations
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
308	Level 1 Isolation Gown, Universal, Yellow	10/bg, 10 bg/cs
308XL	Level 1 Isolation Gown, XL, Yellow	10/bg, 10 bg/cs

Style Guide



Wrist Style



Neck Style



AAMI LEVEL 2

Isolation Gowns

Made of multi-ply SMS material

- Meets AAMI Level 2 requirements for fluid protection
- Proprietary over-the-head design with easy tear away removal
- Elastic wrist with thumb loops to provide a more secure fit under gloves
- Used in low risk situations
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
309	Level 2 Isolation Gown, Universal, Yellow	10/bg, 10 bg/cs
309XL	Level 2 Isolation Gown, XL, Yellow	10/bg, 10 bg/cs

Style Guide



Wrist Style




Neck Style



CoverSoft™ Disposable Apparel

Made of durable and extra soft SMS material that feels like real fabric

Fabric Detail



- Extra soft fabric provides optimal comfort
- Scrub pants have a drawstring tie waist, one front pocket, and one back pocket
- Shirts have two pockets (one top left, one bottom right)
- Dark blue for maximum coverage
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
376S	Shirt, Small	10/bg, 5bg/cs
376M	Shirt, Medium	10/bg, 5bg/cs
376L	Shirt, Large	10/bg, 5bg/cs
376XL	Shirt, X-Large	10/bg, 5bg/cs
376XXL	Shirt, XX-Large	10/bg, 5bg/cs
381S	Pants, Small	10/bg, 5bg/cs
381M	Pants, Medium	10/bg, 5bg/cs
381L	Pants, Large	10/bg, 5bg/cs
381XL	Pants, X-Large	10/bg, 5bg/cs
381XXL	Pants, XX-Large	10/bg, 5bg/cs

Disposable Scrub Wear

Made of durable spunbonded polypropylene



- Scrub pants have a drawstring tie waist, one front pocket, and one back pocket
- Shirts have two pockets (one top left, one bottom right)
- Dark blue
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
375S	Shirt, Small	10/bg, 5bg/cs
375M	Shirt, Medium	10/bg, 5bg/cs
375L	Shirt, Large	10/bg, 5bg/cs
375XL	Shirt, X-Large	10/bg, 5bg/cs
375XXL	Shirt, XX-Large	10/bg, 5bg/cs
380S	Pants, Small	10/bg, 5bg/cs
380M	Pants, Medium	10/bg, 5bg/cs
380L	Pants, Large	10/bg, 5bg/cs
380XL	Pants, X-Large	10/bg, 5bg/cs
380XXL	Pants, XX-Large	10/bg, 5bg/cs

Disposable Coveralls

Made of durable spunbonded polypropylene

Style Guide



Wrist Style

- Elastic wrist and ankle cuffs
- Front zip closure
- Available in two different colors



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
382XL	White, XL	5/bg, 5bg/cs
382XXL	White, XXL	5/bg, 5bg/cs
382XXXL	White, XXXL	5/bg, 5bg/cs
392XL	Blue, XL	5/bg, 5bg/cs



Bouffant Caps

Made of a comfortable spunbonded polypropylene



- Cool and can be worn for hours
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
321B	21", Blue	100/bx, 5 bx/cs
321-10	21", Accordion Style, Blue	100/bg, 10 bg/cs
324B	24", Blue	100/bx, 5 bx/cs
324-10	24", Accordion Style, Blue	100/bg, 10 bg/cs

Shoe Covers

Made of spunbonded polypropylene



- Non-skid to provide better traction
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
350	Shoe Covers, Non-Skid, Blue	100/bx, 3 bx/cs
350-10	Shoe Covers, Non-Skid, Blue	100/bg, 10 bg/cs

Lab Coats

Multi-layer SMS material to provide quality and comfort

- Full length with knit collar and cuff
- Convenient hang tab and size tag sewn into back of coat
- White
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile - Without Pockets		
340	Small	10/bg, 5 bg/cs
341	Medium	10/bg, 5 bg/cs
342	Large	10/bg, 5 bg/cs
343	XL	10/bg, 5 bg/cs
344	XXL	10/bg, 5 bg/cs
345	XXXL	10/bg, 5 bg/cs
Non-Sterile - With Pockets		
340P	Small	10/bg, 5 bg/cs
341P	Medium	10/bg, 5 bg/cs
342P	Large	10/bg, 5 bg/cs
343P	XL	10/bg, 5 bg/cs
344P	XXL	10/bg, 5 bg/cs
345P	XXXL	10/bg, 5 bg/cs



Patient Exam Wear

Made of spunbonded polypropylene

- Made of spunbonded polypropylene
- Gown has open back with belt
- Gown also available in a deluxe version with neck tie
- Shorts have elastic waist and no pockets
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
304	Exam Gown, Dark Blue	10/bg, 5 bg/cs
305	Deluxe Exam Gown, Dark Blue	10/bg, 5 bg/cs
360	Shorts, Unisize, Dark Blue	10/bg, 5 bg/cs



General Use Face Mask

Lightweight mask that provides comfortable face protection



- 3-ply, pleated with earloop
- Not for medical use
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1525	General use, ear loop (not for medical use)	50/bg, 40 bg/cs

ASTM LEVEL 1

Surgical Masks

3-ply pleated lightweight, breathable material



- BFE of $\geq 99\%$
- Meets ASTM Level 1 standards
- Available in earloop or tie on with flexible nosepiece
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1540	Tie, Blue	50/bx, 6 bx/cs
1541	Ear Loop, Blue	50/bx, 6 bx/cs

ASTM LEVEL 2

Surgical Masks

3-ply pleated lightweight, breathable material



- BFE of $\geq 99\%$
- Meets ASTM Level 2 standards
- Available in earloop or tie on with flexible nosepiece
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
CM-101	Ear Loop, Blue	100/bg, 6 bg/cs
CM-102	Tie, Blue	50/bx, 6 bx/cs

Performance Characteristics	Test Standard Methods	ASTM Level 1 Requirements
Bacterial Filtration Efficiency Performance (BFE) (%)	ASTM F2101	$\geq 95\%$
Particulate Filtration Efficiency Performance (PFE) (%)	ASTM F2299	$\geq 95\%$
Differential Pressure (Delta-P) (mm H ₂ O/cm ²)	MIL-M-36945C EN 14683	< 5.0
Fluid Resistance Performance (mmHg)	ASTM 1862	80 mmHg
Flammability	16 CFR part 1610	Class I No Flame Spread

Performance Characteristics	Test Standard Methods	ASTM Level 2 Requirements
Bacterial Filtration Efficiency Performance (BFE) (%)	ASTM F2101	$\geq 98\%$
Particulate Filtration Efficiency Performance (PFE) (%)	ASTM F2299	$\geq 98\%$
Differential Pressure (Delta-P) (mm H ₂ O/cm ²)	MIL-M-36945C EN 14683	< 5.0
Fluid Resistance Performance (mmHg)	ASTM 1862	120 mmHg
Flammability	16 CFR part 1610	Class I No Flame Spread

Protective Plastic Face Shield

Clear polyethylene for optimal visibility

- Double-sided anti-fog
- Protection from splashing or spraying
- Foam headband for gentle comfort
- Single Use
- Not made with natural rubber latex

Item #	Description	Packaging
CM-105	Face Shield	5/bg, 40 bg/cs



Digital Thermometer

Quick 10 Second Read with easy to read large LCD display and low battery indication

- Reads both Centigrade and Fahrenheit; Range: 32° – 43.9°C, 90° – 109.9° F
- Alarm/Beeper when reading is complete
- Battery life: 200 hours of continuous operation
- Auto shut off after 9 minutes of non use
- Built in memory enabling display of the previously measured value

Item #	Description	Packaging
3210D-360	Dual Scale, 10 second read	12 ea/bx, 30 bx/cs



Disposable Stethoscopes

Single Head Stethoscope with 22" Y Piece Tubing

- Overall Length: 32"
- Individually poly bagged
- Single Patient Use
- Not made with natural rubber latex

Item #	Description	Packaging
1115	Yellow	10/bx, 10 bx/cs



Bowles Stethoscope

Chrome plated chest piece with 22" black PVC "Y" tubing and adjustable chrome plated binaural

- Classic blood pressure reading stethoscope
- Lifetime Warranty
- Not made with natural rubber latex

Item #	Description	Packaging
1211	Black	1/bx



Single Head Stethoscopes

Features 22" PVC tubing, anodized aluminum chest piece, and chrome plated brass binaural

- Basic single head stethoscope for blood pressure readings
- Color coded non chill ring
- One year warranty
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1100	Black	1/bx
1100BL	Blue	1/bx
1100GR	Green	1/bx
1100GY	Grey	1/bx
1100R	Red	1/bx



Dual Head Stethoscopes

Features 22" PVC tubing, anodized aluminum chest piece, and chrome plated brass binaural

- Basic dual head stethoscope
- One year warranty
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1200	Black	1/bx
1200BL	Blue	1/bx
1200GR	Green	1/bx
1200GY	Grey	1/bx
1200R	Red	1/bx



Liquid Dispenser

Dispensing design produces steady stream of liquid controlled by finger pressure

- Size: 4-3/4" x 2-1/2"

Item #	Description	Packaging
Non-Sterile		
4026	9 oz. Lotion Style Pump	12 ea/bx, 12 bx/cs




Non-Latex Tourniquets

Used for constricting or compression

- 18"x1" textured nitrile
- Hypoallergenic
- Superior strength
- Convenient packaging designed for kit assembly and resale
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4408	1" x 18", Powder Free	10/bg, 50 bg/bx



Quick Release Tourniquet

Single handed use tourniquet with push button to release

- Universal in size and adjustments will work with child, adult, and thigh size
- Reusable
- Not made of natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4418	Red	12/bx



Latex Finger Cots

Pre-rolled and non-powdered

- Provide user protection and to help prevent contamination and infection
- Approved for medical use and intended to protect health care professionals and patients

Item #	Description	Packaging
Non-Sterile		
4423S	Small	144/bx
4423M	Medium	144/bx
4423L	Large	144/bx
4423XL	X-Large	144/bx



Nitrile Finger Cots

Pre-rolled natural white, powder free, and anti-static

- For medical exam use
- Minimum of .08mm thickness
- Not made of natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4404S	Small	144/bx
4404M	Medium	144/bx
4404L	Large	144/bx
4404XL	X-Large	144/bx



IV Stands

Available in chrome plated steel or stainless steel

- Adjustable from 47" to 84"

Item #	Description	Weight Capacity	Base Diameter	Packaging
4353	Four Caster, 2 Hook	15 pounds	23"	1/ea
4355	Five Caster, 2 Hook	22 pounds	20"	1/ea
4356	Five Caster, 4 Hook	22 pounds	20"	1/ea
4356-SS	Stainless Steel, 4 Hook	22 pounds	20"	1/ea



Body Bag

Made of strong polyvinyl chloride (PVC) fabric

- Bottom of bag is leak proof
- 350lb weight capacity
- Three toe tags



Toe Tags



Item #	Description	Packaging
400B	Body Bag, White	10/cs



Dukal Corporation

2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800

Toll Free: 1 800 243 0741

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Find our full line of products at www.Dukal.com