



Better Health. Better Future.

Dental Care

A commitment to quality solutions that will make you smile.





**Enhancing family &
community health
is the way to a
better future**



Cosmetic Dentistry

Disposables

Endodontics

Evacuation

Impression

Infection Control

Preventative

Cosmetic Dentistry – Micro Applicators & Pre-Bent Needles



Micro Applicator Tips

- Non-linting, non-absorbent fibers suspend material without dripping, thus minimizing waste
- Bendable neck allows access to hard-to-reach areas
- For use with etchants, fluoride varnishing, bonding agents, etc.
- Available in extra slim, super fine, fine, and regular applicator heads
- (4) 100-pc dispenser bottles per box
- Non-sterile
- Single use

Item#	Description	Packaging
UMA-9112-RG	Regular Green	400/BX, 16 BX/CS
UMA-9113-RO	Regular Orange/Peach	400/BX, 16 BX/CS
UMA-9114-RP	Regular Purple	400/BX, 16 BX/CS
UMA-9115-RA	Regular Assorted	400/BX, 16 BX/CS
UMA-9116-FY	Fine Yellow	400/BX, 16 BX/CS
UMA-9117-FP	Fine Pink	400/BX, 16 BX/CS
UMA-9118-FA	Fine Assorted Yellow and Pink	400/BX, 16 BX/CS
UMA-9119-SW	Super Fine White	400/BX, 16 BX/CS
UMA-9014-XS	Extra Slim Black	400/BX, 16 BX/CS

Pre-Bent Dispensing Tips

- Blunt end, pre-bent tips
- Provides fast and accurate dispensing of etchant, resins, sealants, and flowable composite
- Universal Luer Lock design
- Available in color-coded 18G, 20G, 22G, and 25G sizes
- Non-sterile
- Single use



Item#	Description	Packaging
URN-55218P	18G - Pink	100/BX, 20 BX/CS
URN-55220K	20G - Black	100/BX, 20 BX/CS
URN-55220Y	20G - Yellow	100/BX, 20 BX/CS
URN-55222G	22G - Gray	100/BX, 20 BX/CS
URN-55225B	25G - Blue	100/BX, 20 BX/CS

Disposables – Cotton Products



Non-Woven Sponges



- 4-ply non-woven sponge
- High absorbency
- Non-sterile
- Soft, smooth, strong, and virtually non-linting
- Available in premium and economy
- Single use

Item#	Description	Packaging
UNW-4022	2" x 2" Premium Non-Woven	200/SLV, 25 SLV/CS
UNW-3522	2" x 2" Economy Non-Woven	200/SLV, 25 SLV/CS
UNW-3033	3" x 3" Premium Non-Woven	200/SLV, 20 SLV/CS
UNW-3044	4" x 4" Premium Non-Woven	200/SLV, 10 SLV/CS



Cotton-Filled Gauze



- 8-ply 100% white cotton filled
- Soft, smooth, strong & lint-free
- Maximum patient comfort
- Superior absorbency
- No chemical additives
- Non-sterile
- Single use

Item#	Description	Packaging
UCG-1022	2" x 2" Cotton Filled Gauze	200/SLV, 25 SLV/CS
UCG-1044	4" x 4" Cotton Filled Gauze	200/SLV, 10 SLV/CS

Cotton Rolls



- 100% pure white cotton fiber
- High absorbency
- Non-sterile
- No chemical whiteners or additives
- Soft & flexible to increase patient's comfort
- Single use

Item#	Description	Packaging
UCR	Cotton Roll, 1-1/2" x 3/8"	50/RL, 40 RL/BX, 6 BX/CS

Disposable Cotton Tip Applicators



- Made of natural, splinter-free birch wood
- 100% pure white cotton fiber tip
- Each autoclavable bag has a printed indicator for steam process
- Non-sterile
- Single use

Item#	Description	Packaging
UCA-9003	3" Long Wooden Shaft	100/BG, 10 BG/BX, 10 BX/CS
UCA-9006	6" Long Wooden Shaft	100/BG, 10 BG/BX, 10 BX/CS

Disposables – Patient Bibs & Bib Clips



Disposable Patient Bibs

- 2-ply paper / 1-ply poly
- Excellent absorbency
- High quality poly backing keeps patients dry
- Available in a variety of popular colors



Packaging: 500/CS

Available Colors:



Gray
UPB-6099



Blue
UPB-6101



White
UPB-6102



Mauve
UPB-6103



Lavender
UPB-6104



Green
UPB-6105



Yellow
UPB-6106



Peach
UPB-6107



Beige
UPB-6108



Pink
UPB-6109



Black Patient Bibs

- 2-ply paper / 1-ply poly provides excellent absorbency
- Durable, high quality Italian paper
- Vertical embossing and reinforced edge for added strength



Item#

UPB-6110

Description

Black

Packaging

500/CS

Autoclavable Bib / Napkin Holders

- Flexible silicone band provides patient comfort
- Metal clips keep patient bib in place
- Available in 8 popular colors to coordinate with bibs
- Non-sterile



Packaging: 1/BG, 50 BGS/CS

Available Colors:



Light Blue
UBC-6401



White
UBC-6402



Dark Blue
UBC-6403



Lavender
UBC-6404



Green
UBC-6405



Yellow
UBC-6406



Pink
UBC-6409



Red
UBC-6410

Disposables – Tray Covers, Cups & Tongue Depressors



5 oz. Disposable Cups (Plastic and Paper)

- Rolled lip for patient comfort
- Available in a variety of popular colors

Packaging: Plastic: 1000/CS Paper: 800 PCS/CS

Available Colors:



Black



Blue



White



Mauve



Lavender



Green



Yellow



Peach



Beige

PLASTIC:	-	UBC-6201	UBC-6202	UBC-6203	UBC-6204	UBC-6205	UBC-6206	UBC-6207	UBC-6208
PAPER:	UBC-6210	UBC-6211	UBC-6212	-	UBC-6214	UBC-6215	-	-	-



Paper Tray Covers (Ritter)

- Disposable paper cover lays directly on top of tray and assists in preventing cross-contamination
- Specially formulated paper is designed to impede the flow of moisture to the tray below
- 8-1/2" x 12-1/4" with contoured edges perfectly fits standard "B" Trays
- Available in a variety of popular colors

Packaging: 1000/CS

Available Colors:



Blue

UBC-80141



White

UBC-80142



Lavender

UBC-80144



Green

UBC-80145



Yellow

UBC-80146



Pink

UBC-80149




Disposable Tongue Depressors

- Made of natural, splinter-free birch wood with smooth, polished edges
- Highly suitable for oral use
- Non-sterile

Item#	Description	Packaging
UTD-9007	6" x 3/4"	500/BX, 10 BX/CS

New Products – Bendable Brush, Pre-Pasted Toothbrush, Mixing Wells & Dappen Dishes



Bendable Brushes In Slide-Out Box

- Bends to any angle and stays in place
- Innovative slide out box for ease of use
- Handle can be snapped off to shorten
- Can be used for lingual or interproximal use
- Tapered bristles available in fine and regular
- Assorted green, yellow, pink, and blue brushes

Item#	Description	Packaging
UMA-9121-F	Green, Yellow, Pink & Blue Handles w/ Fine White Bristles	100/BX, 16 BX/CS
UMA-9122-R	Green, Yellow, Pink & Blue Handles w/ Regular Black Bristles	100/BX, 16 BX/CS
UMA-9123-A	Assorted Colors w/ Fine White & Regular Black Bristles	100/BX, 16 BX/CS



NEW DISPENSER BOX
Each box contains 100 brushes
Just wet FastBrush to activate and brush as normal

FastBrush™ The Anytime, Anywhere Toothbrush Pre-pasted Toothbrush with Xylitol

- Great for travel, after lunch or before a meeting
- 4 assorted colors
- Individually wrapped
- New dispenser box

Item#	Description	Packaging
UDA-4001	Pre-Pasted Toothbrush	100/BX, 10 BX/CS



Disposable Mixing Wells

- Disposable mixing wells with convenient slot for resting the applicator while not in use
- Economical and lightweight

Item#	Description	Packaging
UMW-9015	Sturdy 2-well Design - A, B Markings	200/BX, 24 BX/CS
UMW-9016	2-well Design - 1, 2 Markings	500/BX, 24 BX/CS
UMW-9017	4-well Design - 1, 2, 3, 4 Markings	200/BX, 5 BX/CS



Disposable Dappen Dishes

- Holds prophylaxis paste, acrylics and amalgams, resins, sealants and other materials
- Tabs bend up to form a handle

Item#	Description	Packaging
UDD-9018-A	Assorted (5) Colors	1000/BX, 10 BX/CS
UDD-90181	Blue	200/BX, 48 BX/CS
UDD-90182	White	200/BX, 48 BX/CS
UDD-90184	Lavender	200/BX, 48 BX/CS
UDD-90186	Yellow	200/BX, 48 BX/CS
UDD-90189	Pink	200/BX, 48 BX/CS

Endodontics – Irrigation Needles & Syringes



Irrigation Syringes (Luer Lock)

- Luer Lock securely holds irrigation syringe tips
- Smooth action plunger for controlled dispensing
- Not for administering injection or dispensing medication
- Available in 3cc, 6cc and 12cc
- Non-sterile
- Single use

Item#	Description	Packaging
URS-55003	3cc	100/BG, 10 BG/CS
URS-55006	6cc	100/BG, 10 BG/CS
URS-55012	12cc	100/BG, 10 BG/CS



Notched Tip Irrigation Needles (Luer Lock)

- 1" bendable "notched" tip permits greater irrigation of tooth canal
- Fluid pressure at apex is significantly reduced by notch tip design
- Luer Lock design allows for easy attachment on syringe
- Fast and easy identification of sizes with color-coded hubs
- Available separately or in a combo pack with syringes
- Comparable to Monoject®
- Non-sterile
- Single use

Item#	Description	Packaging
URN-55223	23G, Blue	100/BX, 20 BX/CS
URN-55227	27G, Yellow	100/BX, 20 BX/CS
URN-55230	30G, Purple	100/BX, 20 BX/CS



Side Port Irrigation Probes (Luer Lock)

- 1" bendable tip with closed blunt end.
- Side-port dispersal creates an upward flushing which allows for easy irrigation of root canals and periodontal pockets.
- Luer Lock design allows for easy attachment on syringe
- Fast and easy identification of sizes with color-coded hubs
- Available separately or in a combo pack with syringes
- Comparable to Max-i-Probe®
- Non-sterile
- Single use


Item#	Description	Packaging
URN-55323	23G, Blue	100/BX, 20 BX/CS
URN-55327	27G, Yellow	100/BX, 20 BX/CS
URN-55330	30G, Purple	100/BX, 20 BX/CS

Endodontics – Needles, Syringes & Organizers



Pre-Tipped Irrigation Syringe Combo Sets (Luer Lock)

- Available with notched tip or side port/probe tip
- Box contains 100 sets of pre-tipped syringes
- Luer Lock securely holds irrigation syringe tips
- Smooth plunger for controlled dispensing
- Not for administering injection or dispensing medication
- Syringe holds 3cc
- Non-sterile
- Single use



NOTCHED NEEDLE TIP: Monoject® Style

Item#	Description	Packaging
URN-55003-23	23 Gauge Blue Hub/Notched	100 SETS/BX, 20 BX/CS
URN-55003-27	27 Gauge Yellow Hub/Notched	100 SETS/BX, 20 BX/CS
URN-55003-30	30 Gauge Purple Hub/Notched	100 SETS/BX, 20 BX/CS

SIDE PORT / PROBE NEEDLE TIP: Max-i-Probe® Style

Item#	Description	Packaging
URN-55103-23	23 Gauge Blue Hub/Side Port	100 SETS/BX, 20 BX/CS
URN-55103-27	27 Gauge Yellow Hub/Side Port	100 SETS/BX, 20 BX/CS
URN-55103-30	30 Gauge Purple Hub/Side Port	100 SETS/BX, 20 BX/CS

Endo Organizer Rings

- Includes 1 organizer and 2 rings (Large and Small)
- Kit also comes with 50pcs of disposable foam refills
- Built-in ruler for easy measurement of files up to 30mm
- Can be worn on either right or left hand
- Ring and holder are easily cleaned
- Non-sterile
- Single use



Item#	Description	Packaging
URE-10200G	2 Rings + 50 Green foam inserts	20 BG/CS

Evacuation – HVE, Aspirator Tips, Canisters & Saliva Ejectors







High Volume Evacuation Tips

- Designed to fit all high-volume evacuation systems
- Vented on one side / non-vented on the other
- Made of rigid opaque plastic
- Available in a variety of popular colors
- Non-sterile
- Single use

Packaging: 100/BG, 10 BG/CS

Available Colors:



Surgical Aspirator Tips

- Designed to fit all high-volume evacuation systems
- 30° angle for maximum patient comfort
- Made of rigid opaque plastic
- Available in 1/4", 1/8", and 1/16"
- Non-sterile
- Single use

Item#	Description	Packaging
UST-9073	1/8" x 6-1/4" (White)	25/BG, 10 BG/CS
UST-9074	1/4" x 6-1/4" (Green)	25/BG, 10 BG/CS
UST-9075	1/16" x 7-1/4" (Blue)	25/BG, 10 BG/CS



UNiVAQ Disposable Canisters

- Disposable canister, lid, mesh filter and gasket are designed to replace the present canister and filter in your central vacuum pump.
- Clear plastic for easy viewing
- Saves staff time on cleaning
- Minimizes contact with infectious material
- Non-sterile




Item#	Description	Packaging
UBC-82200	Evac Canisters with Mesh Screen 2-3/4" x 3-5/8"	12/BX, 6 BX/CS
UBC-82300	Evac Canisters with Mesh Screen 3-1/2" x 4-3/8"	8/BX, 6 BX/CS
UBC-82350	Evac Canister with Finer Mesh Screen 3.5"x4-3/8"	8/BX, 6 BX/CS

Evacuation – HVE, Aspirator Tips, Canisters & Saliva Ejectors



UBC-85500 (5500) Lavender
Fits: A-dec, Cascade, and Performer
Diameter: 1-1/4"




UBC-85501 (5501) Blue
Fits: A-dec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus, MDT, Schein **Diameter: 2-1/8"**



UBC-85502 (5502)
Fits: A-dec models with metal canisters
Diameter: 2-1/4"




UBC-85503 (5503) Green
Fits: Ampco, Belmont, Biotec, Chayes, Dabi-Atlante, Dansereau, DCI, Den-Tal-Ez, Dentech, Forest, Knight, Pelton & Crane, Proma
Diameter 1-7/8"



UBC-85505 (5505)
Fits: Pelton & Crane Spirit with white plastic canisters
Size: 2' x 2-3/4"



UBC-85508S (5508-S)
Fits: All saliva ejectors
Universal fit
*72 pcs/bag - 50 bags/case



UBC-85509 (5509)
Fits: DCI/Belmont Reality, Series, System 6000.
Diameter: 1-1/2"



Disposable Evacuation Traps

- Directly replaces current evacuation or cuspidor trap
- Disposable and easily replaceable
- Made of rigid white plastic, except where noted
- Available in a variety of sizes to fit your needs

Packaging: 144/BX, 8 BX/CS *Except where noted



UBC-8HVE (HVE-CAP)
Fits: On all standard high volume evacuator tips
Universal fit
* 144/BG - 15 BG/CS



UBC-85512 (5512)
Fits: A-dec 500 series
Diameter: 2-3/4"



UBC-86100C (6100-C)
Fits: Most cuspidors, including A-dec, Ampco, Belmont, Dabi-Atlante, Dansereau, DCI, Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter, Schein
Universal fit



UBC-86200C (6200-C)
Fits: Siemens, Planmeca, and other cuspidors with small openings
Universal fit



UBC-86400-C
Fits: A-dec Cascade, Performer
Diameter: 1-5/8"

Impression – Alginate

CAVEX ColorChange Fast Set Alginate



- Indicates three color changes
- 9 days of dimensional stability
- Tear and deformation-resistant
- Detail reproduction of at least 25 µm
- Maximum compatibility with plaster
- Short setting time
- Highly elastic
- Dust-free



MIXING SCHEME:

Mixing	Filling+Placing	Setting in Mouth
30 sec.	1 minute	1 minute
Purple	Light Pink	White



CLINICAL RATING OF 96%
BY THE DENTAL ADVISOR

Item#	Description	Packaging
AA323	Color Change, Fast Set	500G/PK, 20 PK/CS

CAVEX Normal Set & Fast Set Alginate




- Solid alginate impression
- Detail reproduction of 50 µm
- Optimal tear strength and deformation resistant
- Natural pink color
- Peppermint flavor
- Detail reproduction of at least 50 µm
- Compatible with all plaster types
- Short setting time
- Dust-free
- Available in Normal and Fast Set

Item#	Description	Packaging
AA025	CA37 Fast Set (setting time in mouth 1 min.)	500G/PK, 20 PK/CS
AA075	CA37 Normal Set (setting time in mouth 1-1/2 min.)	500G/PK, 20 PK/CS

AA075
CA37 Normal Set

AA025
CA37 Fast Set

CAVEX Orthotrace Alginate



- Clearly legible pink color
- Pleasant fruit aroma
- Extra fast setting
- Rapid water absorption
- Highly elastic
- Short Setting time
- Dust-free

Item#	Description	Packaging
AA380	Orthodontrace Extra Fast Set	500G/PK, 20 PK/CS

**extra fast set (2.10 min incl. mixing) for efficient and fast working

CAVEX Cream Normal & Fast Set



- Detail reproduction of at least 5 µm
- Easy mixing
- Highly elastic
- High tear strength
- Five-year shelf life
- Outstanding scan ability
- Dust-free
- Available in Normal and Fast Set

Item#	Description	Packaging
AA605	Cream Fast Set	500G/PK, 20 PK/CS
AA600	Cream Normal Set	500G/PK, 20 PK/CS



CLINICAL RATING OF 97%
BY THE DENTAL ADVISOR

Impression – Alginate

Temporary Cement

- Easy to use
- Eugenol free
- Distinctive color
- Perfect sealing
- Low solubility
- Short curing time
- Easy to remove



Item#	Description
CE090	1 Tube 35g Base, 1 Tube 16g Catalyst and 1 Mixing Pad



#AT021



#AT076

#AT042



#AT041



#AT240


Cavex Alginate Accessories

Item#	Description	Packaging
AT021	Water Dosing Bottle	1 BTL/BX
AT041	Flexible Mixing Bowl	1 PC/EA
AT042	Flat Mixing Spatula	1 PC/EA
AT240	Blue Storage Container	1 SET/BX
AT260	Pink Storage Container	1 SET/BX
AT076	Scoop and Water Measuring Cup (Two sets provided for each case of Alginate purchased)	2 SET/CS

Impression – Impression Accessories

Cavex alginates have long exceeded the criteria established by institutions such as ISO/ADA for this process. Cavex has set its own standard, which is reflected in the Cavex alginate branding.

Cavex Alginates, best in test



Tensile strength of rival alginate brands versus Cavex alginate

Sixty years of research and experience have resulted in superior alginate impression materials with an unblemished worldwide reputation.

Every day, tens of thousands of dental professionals benefit from the unparalleled quality and advantages of the Cavex alginate system.

Exactly why are Cavex alginates so good that they leave the competition far behind?

It's because impression materials, and alginates in particular, are Cavex's core business and so Cavex strives for absolute perfection in this area.

Cavex alginates – performance and specifications

	Cavex CA37	Cavex ColorChange	Cavex Orthotrace	Cavex Cream Alginate
Contents	500g.	500g.	500g.	500g.
Number of imprints	approx. 35	approx. 38	approx. 35	approx. 35
Consistency	solid	elastic	solid	elastic
Dimensional stability	2 days	9 days	2 days	2 days
Detail reproduction	50µm	25µm	50µm	5µm
Colour	pink	light pink	magenta	purple
Colour change	no	yes	no	no
Fragrance	peppermint	peppermint	red fruit	bubble gum
Setting types	NS & FS	FS	Extra FS	NS & FS
Shelf life	5 years	5 years	5 years	5 years
Mixing time in sec.	NS = 30 / FS = 30	30	30	NS = 30 / FS = 20
Curing time in mouth in sec.	NS = 90 / FS = 60	60	60	NS = 90 / FS = 60
Total curing time in sec.	NS = 210 / FS = 150	150	130	NS = 210 / FS = 140



Tear-Resistant and Strong

- 2.5 x ISO/ADA criterion as a minimum
- Can be cast at least twice
- No re-working



Snap Set

- Predictable processing time
- Fully cured in a short time
- Minimal patient stress



5-year shelf life

- Developed to perfection
- Highest quality raw materials
- Outstanding stability
- Perfect impressions for five years

Impression – Bite Trays



The tray can be heated and molded for a more custom fit.



Tray hole openings can be easily extended with a dental bur to accept implants.



COLOR-Brite

Perforated Transparent Impression Trays

- Disposable premium impression trays made of high quality, rigid, transparent plastic
- Strategically perforated trays allow easy modification to accommodate for implants
- Trays can be heated and molded for a custom fit
- Numbered handles for easy position identification
- Easy to apply I.D. labels included in every bag
- Non-sterile
- Single use

Item#	Description	Packaging
UTA-7201	#1 Large: Upper, Yellow	12/BG, 10 BG/CS
UTA-7202	#2 Large: Lower, Yellow	12/BG, 10 BG/CS
UTA-7203	#3 Medium: Upper, Blue	12/BG, 10 BG/CS
UTA-7204	#3 Medium: Lower, Blue	12/BG, 10 BG/CS
UTA-7205	#5 Small: Upper, Orange	12/BG, 10 BG/CS
UTA-7206	#6 Small: Lower, Orange	12/BG, 10 BG/CS
UTA-7207	#7 Partial: Upper Left/Lower Right, Green	12/BG, 10 BG/CS
UTA-7208	#8 Partial: Upper Right/Lower Left, Green	12/BG, 10 BG/CS
UTA-7209	#9 Anterior: Upper, Pink	12/BG, 10 BG/CS
UTA-7210	#10 Anterior: Lower, Pink	12/BG, 10 BG/CS
UTA-7210-A	Assorted, 4 Pcs Medium & 2 Pcs of each other size	12/BG, 10 BG/CS

BLU-Tray

Perforated Impression Trays

- Made of high quality, rigid plastic
- Perforated disposable impression trays allow a perfect impression at an economical cost
- Non-sterile
- Single use



Item#	Description	Packaging
UTA-7211	#1 Large: Upper, Yellow	12/BG, 10 BG/CS
UTA-7212	#2 Large: Lower, Yellow	12/BG, 10 BG/CS
UTA-7213	#3 Medium: Upper, Blue	12/BG, 10 BG/CS
UTA-7214	#3 Medium: Lower, Blue	12/BG, 10 BG/CS
UTA-7215	#5 Small: Upper, Orange	12/BG, 10 BG/CS
UTA-7216	#6 Small: Lower, Orange	12/BG, 10 BG/CS
UTA-7217	#7 Partial: Upper Left/Lower Right, Green	12/BG, 10 BG/CS
UTA-7218	#8 Partial: Upper Right/Lower Left, Green	12/BG, 10 BG/CS
UTA-7219	#9 Anterior: Upper, Pink	12/BG, 10 BG/CS
UTA-7220	#10 Anterior: Lower, Pink	12/BG, 10 BG/CS
UTA-7220-A	Assorted, 4 Pcs Medium & 2 Pcs of each other size	12/BG, 10 BG/CS

Impression – Bite Trays & Syringes



3-Way Impression Bite Trays



- Disposable premium bite trays precisely capture the master, counter impression and bite registration in a single procedure
- U-HOLD grooves provide superior material retention
- Rigid frame with nylon webbing to ensure accurate bite
- Ergonomically designed thumb grip aids in placement and removal
- Non-sterile
- Single use

Item#	Description	Packaging
UBT-7311	Full Arch - Pink	30/BX, 10 BX/CS
UBT-7312	Anterior- Lavender	35/BX, 10 BX/CS
UBT-7313	Quadrant- Yellow	35/BX, 10 BX/CS
UBT-7314	Posterior-Blue	50/BX, 10 BX/CS
UBT-7315	Sideless Posterior- White	50/BX, 10 BX/CS
UBT-7316	Assorted, 7 Pcs for Each Size	35/BX, 10 BX/CS




Curved Tip Utility Syringes



- Ideal for pressure irrigation and impression material injection
- 12cc
- Non-sterile
- Single use

Item#	Description	Packaging
URS-55112	12cc	50/BX, 10 BX/CS

Disposable Impression Syringes



- Tip is bendable
- For use with all types of impression materials
- Seamless construction prevents splitting during extrusion
- Smooth sliding plunger
- Non-sterile
- Single use

Item#	Description	Packaging
URS-55900	5ml	50/BG, 10 BG/CS


Oral Tips



- Core oral tips for use with build-up material
- Intra oral tips attach securely and provide a precise dispensing of the impression material
- Non-sterile
- Single use

Item#	Description	Packaging
UMT-9022	Intra Oral Tips, White	100/BG, 10 BG/CS
UMT-9023	Core Oral Tips, Short	25/BG, 10 BG/CS
UMT-9024	Core Oral Tips, Long	25/BG, 10 BG/CS


See-Through and Self-Sealing Sterilization Pouches



Item#	Description	Common Cross Reference Size	Packaging
UPP-223	2-1/4" x 3-3/4"	(2-1/4" x 2-3/4")	200/BX, 30 BX/CS
UPP-225	2-1/4" x 5"	(2-1/4" x 4")	200/BX, 16 BX/CS
UPP-2710	2-3/4" x 10"	(2-3/4" x 9")	200/BX, 16 BX/CS
UPP-2710-1K	2-3/4" x 10"	(2-3/4" x 9")	1000/BX, 4 BX/CS
UPP-29	2" x 9"	(2" x 8")	200/BX, 16 BX/CS
UPP-3565	3-1/2" x 6-1/2"	(3-1/2" x 5-1/4")	200/BX, 16 BX/CS
UPP-3510	3-1/2" x 10"	(3-1/2" x 9")	200/BX, 16 BX/CS
UPP-3510-1K	3-1/2" x 10"	(3-1/2" x 9")	1000/BX, 4 BX/CS
UPP-3212	3-1/4" x 12"	(3-1/4" x 11")	200/BX, 16 BX/CS
UPP-4212	4-1/4" x 12"	(4-1/4" x 11")	200/BX, 12 BX/CS
UPP-5275	5-1/4" x 7-1/2"	(5-1/4" x 6-1/2")	200/BX, 16 BX/CS
UPP-5211	5-1/4" x 11"	(5-1/4" x 10")	200/BX, 12 BX/CS
UPP-7513	7-1/2" x 13"		200/BX, 10 BX/CS
UPP-1015	10" x 15"		200/BX, 8 BX/CS
UPP-1218	12" x 18"		200/BX, 5 BX/CS

- A Dual External Indicators**
Easy to read "Color Changing Indicators" printed on each pouch assuring proper sterilization process by Steam and by Ethylene Oxide.
- B Built-in Internal Indicator**
For Steam Process – to ensure the proper temperature is reached inside the pouch.
- C Complete Film/Paper Separation**
Clean separation (peeling) of the 4-ply see-thru film from the medical grade paper maintains sterility of contents until their use.
- D Blue-Tinted Transparent Film**
For easier identification of contents and possible punctures by contrast of color.
- E Triple Seal & Heat Processed**
For a seal that is consistent and strong which eliminates shattering and tearing.
- F Wide Self-Seal Adhesive Strip**
Easy-seal flap to ensure the integrity of contents after sterilization process.

EASY TO USE



1 Place item in pouch



2 Peel off the paper strip to expose adhesive band



3 Fold over flap and press firmly



Use them once,
Use them forever.

ASA Saliva Ejectors

- Vented and smooth tip for maximum suction and patient comfort
- Soft, non-removable tip for patient safety
- Wire reinforced tube maintains desired shape
- Available in today's most popular colors or with a "fresh Mint" scented body



Item#	Description	Packaging
2900E-BBS15	White/White	100/BG, 10 BG/CS
2900E-TCS15	Clear/Light Blue	100/BG, 10 BG/CS
2900E-PPS15	Pink/Pink	100/BG, 10 BG/CS
2900E-ATS15	Transparent Royal Blue/Clear	100/BG, 10 BG/CS
2900E-ZTS15	Transparent Aqua Blue/Clear	100/BG, 10 BG/CS
2900E-LTS15	Transparent Lavender/Clear	100/BG, 10 BG/CS
2900E-VTS15	Transparent Green/Clear	100/BG, 10 BG/CS




Image Coming Soon

Cut to Length Sterilization Tubing

- Same quality as our Sterilization Pouch, but with the added convenience of being able to be cut to your desired length.
- Easy to read “Color Changing Arrows” printed along the paper side, assures package has been processed by Steam or EO.
- Blue Tinted Transparent Film for easier identification of contents and possible punctures by contrast of color.

Item#	Description	Packaging
UPT-002	2" x 100'	1 RL/BX, 10 BX/CS
UPT-003	3" x 100'	1 RL/BX, 10 BX/CS
UPT-004	4" x 100'	1 RL/BX, 10 BX/CS
UPT-006	6" x 100'	1 RL/BX, 10 BX/CS



When processed, the indicator stripes change color from green to black



Autoclave Indicator Tape

- Printed indicator stripes on tape change color when exposed to the steam sterilization process.
- Specially formulated adhesive adheres to and seals all types of wraps, tubing, pouches, etc. – without the worry of it sliding or “popping off”.
- Pliable paper surface allows you to write on tape with an indelible pen for easy identification and maintaining better records.

Item#	Description	Packaging
UPST-050	½" x 100'	12 RL/BX, 6 BX/CS
UPST-075	¾" x 100'	8 RL/BX, 6 BX/CS
UPST-100	1" x 100'	6 RL/BX, 6 BX/CS

Sterilization CSR Wrap

- Soft, manageable, no memory wrap with high tensile strength.
- Consistent blue texture for easy identification of punctures or tears.
- For use with Steam or EO sterilization processes.
- Repels water, blood, iodine and saline.



Item#	Description	Packaging
UPP-1212	12" x 12"	100/BG, 10 BG/CS
UPP-1515	15" x 15"	100/BG, 10 BG/CS
upp-2020	20" x 20"	100/BG, 10 BG/CS
upp-2424	24" x 24"	100/BG, 10 BG/CS

Infection Control – Disposable Equipment Covers

Easy Peel Adhesive Covers

- 4" x 6" perforated sheets with non-adhesive side edges for easy removal
- Specially formulated adhesive leaves no residue behind
- Helps to prevent cross-contamination
- Convenient dispenser box contains 1 roll of 1200 perforated 4" x 6" sheets



*UBC-8045

Item#	Description	Packaging
UBC-8045*	4" x 6", 1200 sheets/roll - Black	1 RL/BX, 8 BX/CS
UBC-8046	4" x 6", 1200 sheets/roll - Green	1 RL/BX, 8 BX/CS
UBC-8047	4" x 6", 1200 sheets/roll - Pink	1 RL/BX, 8 BX/CS
UBC-8048	4" x 6", 1200 sheets/roll - Blue	1 RL/BX, 8 BX/CS
UBC-8049	4" x 6", 1200 sheets/roll - Clear	1 RL/BX, 8 BX/CS

*Sold under the BLAKCAT brand name

Available Colors:



Chair & Headrest Covers

- Disposable clear poly protects operatory chair from cross-contamination
- Made of a resilient, yet supple clear poly material that resists tearing
- Assists in protecting chair and prolonging the life of your investment
- Available in a variety of sizes for a perfect fit every time



Image Coming Soon

Item#	Description	Packaging
UBC-8021	29" x 80", Full Chair Cover	125/BX, 3 BX/CS
UBC-8022	48" x 56", Full Chair Cover	150/BX, 3 BX/CS
UBC-8023	27 1/2" x 24", Half Chair Cover	225/BX, 3 BX/CS
UBC-8123-R	27 1/2" x 24", Half Chair Cover on a Roll	225/BX, 3 BX/CS
UBC-8024	11" x 9 1/2 x 2" Headrest Cover	250/BX, 12 BX/CS
UBC-8025	14" x 9 1/2 x 2" Headrest Cover	250/BX, 12 BX/CS

Tissue Headrest Covers

- 2-Ply Tissue + Poly
- Sizes 10" x 10" and 10" x 13"
- Available in a variety of popular colors
- Non-sterile
- Single use

Packaging: 500/CS

Available Colors:



Blue



White



Lavender



Green



Pink

Size:

10" x 10"

10" x 13"

UBC-80261

UBC-80271

UBC-80262

UBC-80272

UBC-80264

UBC-80274

UBC-80265

UBC-80275

UBC-80269

UBC-80279

Infection Control – Barrier Covers



Equipment Barriers

- Made of a clear poly material that protects instruments from cross-contamination
- Minimizes staff time spent on cleaning equipment
- Specially designed barriers effectively cover equipment with a near custom fit
- Syringe and pen sleeves have pre-cut openings for easy application

Item#	Description	Packaging
UBC-8028	2-1/2" x 10", Syringe Sleeve w/ Opening, Clear	500/BX, 36 BX/CS
UBC-8029	2-1/2" x 10", Syringe Sleeve w/ Opening, Blue	500/BX, 36 BX/CS
UBC-8030	2-1/2" x 10", Digital X-Ray Sensor Protector	500/BX, 36 BX/CS
UBC-8031	4" x 5-3/4", T-Handle Cover (T-Shape)	500/BX, 60 BX/CS
UBC-8032	23" x 31", Universal X-Ray Cover	200/BX, 3 BX/CS
UBC-8033	15" x 26", Universal X-Ray Cover	250/BX, 6 BX/CS
UBC-8034	Curing Light, Full Pistol	500/BX, 10 BX/CS
UBC-80341-F	2" x 12-1/2" Curing Light, Pen (Small)	500/BX, 10 BX/CS
UBC-80342-S	2" x 12-1/2" Curing Light, Pen (Small)	500/BX, 10 BX/CS
UBC-80343-L	3-1/8" x 12-1/2" Curing Light, Pen (Large)	500/BX, 10 BX/CS
UBC-8035	1-1/2" x 12 1/2" Low Speed Pen Sleeve	500/BX, 36 BX/CS
UBC-8036	1" x 9" High Speed Universal	500/BX, 36 BX/CS
UBC-8037	0 (Small) Sensor Sleeve Universal	500/BX, 36 BX/CS
UBC-8038	2 (Large) Sensor Sleeve Universal	500/BX, 36 BX/CS
UBC-80392	2" x 1200' Long Tube Sleeve	1 RL/BX, 4 RL/CS
UBC-80394	4" x 1200' Long Tube Sleeve	1 RL/BX, 4 RL/CS

Phosphor Plate Barrier Envelopes

SLDR EASY LOAD AND SEAL BARRIER ENVELOPES

SLDR tab allows quick and easy loading and sealing. Once peeled off, the envelope is easily sealed.


- "V" notch for easy opening
- Soft, pliable vinyl with rounded corners
- Transparent side helps in identifying plate location
- Non-sterile • Single use
- Designed to fit all Dentsply, Kodak, Sirona and Air Techniques phosphor plates.



Item#	Description	Size	Packaging
UBE-8160	V-Notch, Top Opening	0	100/BX, 36 BX/CS
UBE-8161	V-Notch, Top Opening	1	100/BX, 36 BX/CS
UBE-8162	V-Notch, Top Opening	2	300/BX, 8 BX/CS

TOP OPENING - BARRIER ENVELOPES

- "V" notch for quick and easy opening
- Soft, pliable vinyl with rounded corners
- Transparent side helps in identifying plate location
- Non-sterile • Single use
- Designed to fit all Dentsply, Kodak, Sirona and Air Techniques phosphor plates.



Item#	Description	Size	Packaging
UBE-8050	V-Notch, Top Opening	0	100/BX, 36 BX/CS
UBE-8051	V-Notch, Top Opening	1	100/BX, 36 BX/CS
UBE-8052	V-Notch, Top Opening	2	100/BX, 36 BX/CS
UBE-8053	V-Notch, Top Opening	3	100/BX, 36 BX/CS

Infection Control – Barrier Covers & Napkin Holders



Tray Sleeves

- Back-flip closure ensures tray is completely covered
- Made of a clear poly material that protects tray from cross-contamination
- Minimizes time spent on cleaning trays
- Available in tray sizes A, B, E, and F

Item#	Description	Size	Packaging
UBC-8011A	A Tray Sleeve	11-5/8" x 14-1/2"	500/BX, 6 BX/CS
UBC-8012B	B Tray Sleeve	10-1/2" x 14"	500/BX, 6 BX/CS
UBC-8010E	E Tray Sleeve	11-1/2" x 16"	500/BX, 6 BX/CS
UBC-8013F	F Tray Sleeve	7-1/2" x 10-1/2"	500/BX, 12 BX/CS




U-Hold Elasticized Bib Holders

- Compact box allows easy storage and economical shipping
- Made of soft, silk-like material for patient comfort
- Can be stretched to fit all patients
- Easy-stick, self-adhesive tabs
- Tangle-free dispenser box
- Reduces the risk of cross-contamination
- More cost-effective than autoclave bib holders
- Non-sterile
- Single use



Item#	Description	Packaging
UBC-6311	Elastic Bib Holder - Blue	250/BX, 30 BX/CS
UBC-6312	Elastic Bib Holder - White	250/BX, 30 BX/CS

NEATER*



SOFTER*



STRETCHIER*



CHEAPER*

*Than the leading brand

Infection Control – Apparel











Lab Coats & Jackets

- Unique PATENT-PENDING 3-position adjustable waist band design provides a truly custom fit
- Knit cuffs and collar
- Snap front
- 2 hip & 1 chest pocket
- Premium SMS material, soft, breathable & fluid resistant
- Washable and autoclavable
- Not made with natural rubber latex

Packaging: 10/BG – 5 BG/CS











Lab Jackets • Same features as lab coats but 25% shorter length

								
Size	Black*	Sky Blue	White	Ciel Blue	Purple/ Lavender	Teal/Green	Raspberry Pink	Bubblegum Pink
Small	UGJ-6500-S	UGJ-6501-S	UGJ-6502-S	UGJ-6503-S	UGJ-6504-S	UGJ-6505-S	UGJ-6509-S	UGJ-6510-S
Medium	UGJ-6500-M	UGJ-6501-M	UGJ-6502-M	UGJ-6503-M	UGJ-6504-M	UGJ-6505-M	UGJ-6509-M	UGJ-6510-M
Large	UGJ-6500-L	UGJ-6501-L	UGJ-6502-L	UGJ-6503-L	UGJ-6504-L	UGJ-6505-L	UGJ-6509-L	UGJ-6510-L
X-Large	UGJ-6500-XL	UGJ-6501-XL	UGJ-6502-XL	UGJ-6503-XL	UGJ-6504-XL	UGJ-6505-XL	UGJ-6509-XL	UGJ-6510-XL
2X-Large	UGJ-6500-XXL	UGJ-6501-XXL	UGJ-6502-XXL	UGJ-6503-XXL	UGJ-6504-XXL	UGJ-6505-XXL	UGJ-6509-XXL	UGJ-6510-XXL



Lab Coats • Same features as lab jackets but 25% longer length

								
Size	Black*	Sky Blue	White	Ciel Blue	Purple/ Lavender	Teal/Green	Raspberry Pink	Bubblegum Pink
Small	UGC-6600-S	UGC-6601-S	UGC-6602-S	UGC-6603-S	UGC-6604-S	UGC-6605-S	UGC-6609-S	UGC-6610-S
Medium	UGC-6600-M	UGC-6601-M	UGC-6602-M	UGC-6603-M	UGC-6604-M	UGC-6605-M	UGC-6609-M	UGC-6610-M
Large	UGC-6600-L	UGC-6601-L	UGC-6602-L	UGC-6603-L	UGC-6604-L	UGC-6605-L	UGC-6609-L	UGC-6610-L
X-Large	UGC-6600-XL	UGC-6601-XL	UGC-6602-XL	UGC-6603-XL	UGC-6604-XL	UGC-6605-XL	UGC-6609-XL	UGC-6610-XL
2X-Large	UGC-6600-XXL	UGC-6601-XXL	UGC-6602-XXL	UGC-6603-XXL	UGC-6604-XXL	UGC-6605-XXL	UGC-6609-XXL	UGC-6610-XXL

Infection Control – Apparel



Isolation Gowns

- Made of high quality and heavyweight SMS material
- Knit cuffs with tie back
- Not made with natural rubber latex
- Non-sterile • Single use

Packaging: 10/BG, 5 BG/CS

Size	Blue	Teal/Green	Yellow	Pink	White
Small	UGI-6711-S	–	–	–	–
Regular	UGI-6711	UGI-6705*	UGI-6707*	UGI-6719	UGI-6712
X-Large	UGI-6711-XL	UGI-6705-XL*	–	–	UGI-6712-XL

*Made of PP material



Anti-Fog Full Face Shields

- Soft foam front allows use of mask and eyeglasses
- No-fog coating provides a crystal clear view
- Comfortable elastic band holds face shield in place
- One size fits all
- Non-sterile
- Single use

Item#	Description	Packaging
ULM-6091	Full Face	200/CS



Non-Woven Caps

- Made of high quality non-woven material
- Sewn-in elastic band holds cap snug and in place
- Available in Blue and White
- Non-sterile
- Single use

Item#	Description	Packaging
UGA-6802	21" Blue	100/BG, 20 BG/CS
UGA-6803	21" White	100/BG, 20 BG/CS
UGA-6804	24" Blue	100/BG, 20 BG/CS
UGA-6805	24" White	100/BG, 20 BG/CS



Non-Woven Shoe Covers

- Made of high quality non-woven material
- Anti-skid bottoms
- Sewn-in elastic band holds covers snug and in place
- Available in Ciel Blue
- Non-sterile
- Single use

Item#	Description	Packaging
UGA-6801	Shoe Cover, Ciel Blue	100/BG, 10 BG/CS

Infection Control – Ear-Loop Face Masks



ASTM F2100
LEVEL 1 2 3

ASTM RATED LEVEL 1 AND LEVEL 3 Black Ear-Loop Face Masks

Level 1

- Meets and/or exceeds ASTM Level 1 Standards
- 4-Layer, 3-fold
- Superior breathability
- Ultra soft inner layer



ASTM F2100
LEVEL 1 2 3

Level 3

- Meets and/or exceeds ASTM Level 3 Standards
- 4-Layer, 4-fold
- Superior breathability
- Ultra soft inner layer

Item#	ASTM Rated	Packaging
ULM-6180	Level 1	50/BX, 20 BX/CS
ULM-6380	Level 3	50/BX, 20 BX/CS

Flex Mask Level 1

- Meets and/or exceeds ASTM Level 1 Standards
- 3-Layer, 4-fold
- Superior breathability
- Ultra soft inner layer
- Pop-out flex center strip

Packaging: 50/BX, 20 BX/CS



a unique **pop-out center strip** to enhance breathability and reduce mask movement when speaking



Blue

ULM-6181



White

ULM-6182



Lavender

ULM-6184



Green

ULM-6185



Yellow

ULM-6186



Pink

ULM-6189

FluidBloq Mask Level 3

- Meets and/or exceeds ASTM Level 3 Standards
- 4-Layer, 4-fold
- Superior breathability
- Ultra soft inner layer
- FluidBloq mesh layer provides superior fluid resistance

Packaging: 50/BX, 20 BX/CS



Blue

ULM-6381



White

ULM-6382



Pink

ULM-6389

Infection Control – Air / Water Syringe Tips




COLOR-Brite

Disposable Air/Water Syringe Tips

- Beveled edge for maximum patient comfort
- Two grooves, one for locking indication, the other for securing a tight lock
- Convenient re-sealable pouch contains an assorted arrangement of colors
- Quad core air-channel for ultimate performance
- Rounded to minimize wear and tear off O-ring in hand piece
- No adapter required
- Non-sterile
- Single use

Item#	Description	Packaging
UAT-9053	Assorted Color Tips	250/BG, 16 BG/CS



Disposable Air/Water Syringe Tips

- Rounded edge for maximum patient comfort
- Convenient dispenser bottle keeps work area clear and assists in preventing cross-contamination
- Every box contains 2 dispenser bottles with 125 pieces each
- Rounded to minimize wear and tear of O-ring in hand piece
- Body made of rigid white plastic with either a white plastic or metal core
- No adapter required
- Non-sterile
- Single use

Item#	Description	Packaging
UAT-9051	Plastic with Metal Core	250/BG, 16 BG/CS
UAT-9052	Plastic with Plastic Core	250/BG, 16 BG/CS



Retainer Box / Denture Box

- One piece clamshell design with easy to open/close “snap-lock”
- Sturdy, colorful plastic helps protect retainer/dentures from dust and damage
- Vented top and bottom holes allow box to “breathe”
- Flat bottom provides area for writing “Return if found” information
- Denture box has deep soaking well to hold both uppers and lowers
- Assorted colors

Item#	Description	Packaging
UDB-7405	Retainer Box	12/BG, 20 BG/CS
UDB-7406	Denture Box	12/BG, 20 BG/CS

Intra-Oral Camera Covers and X-Ray Sensor Sheaths



- Soft and pliable
- Fits tightly over the lens
- Reduces the need for disinfectants
- Custom-fitted with paper backing

Item#	Description	Size	Packaging
UBC-820831	Intra-Oral Camera Cover	Pro-Den Systems/Dent-X, Pro-scope 1000 & Oral Scan, Easy Doc	100/BX, 15 BX/CS
UBC-820855	Intra-Oral Camera Covers	Siemens, Ceracam/Minicam Ultra	100/BX, 15 BX/CS
UBC-820963	Intra-Oral Camera Covers	Video Dental Concepts, Quickcam Smile	100/BX, 15 BX/CS
UBC-821013	Intra-Oral Camera Covers	Digital Doc, Iris	100/BX, 15 BX/CS
UBC-820824	X-Ray Sensor Sheath	Schick, Size 1	100/BX, 15 BX/CS
UBC-820825	X-Ray Sensor Sheath	Schick, Size 2	100/BX, 15 BX/CS
UBC-820861	X-Ray Sensor Sheath	Regam, Size 2	100/BX, 15 BX/CS
UBC-820978	X-Ray Sensor Sheath	Carestream/Kodak 6100 Size 1	100/BX, 15 BX/CS
UBC-820979	X-Ray Sensor Sheath	Carestream/Kodak 6100 Size 2	100/BX, 15 BX/CS
UBC-820999	X-Ray Sensor Sheath	Dexis/Universal	100/BX, 15 BX/CS

New Products – Piksters Interdental Brushes and Protective Eyewear Interdental Brushes



Interdental Brushes with Bamboo Handles

- Bamboo brush handles create less plastic waste
- 100% recycled biodegradable packaging
- Sustainably sourced handles are up to 97% biodegradable

Item#	Description	Packaging
PKB-8-00	Size 00	8/PK, 10 PK/BX
PKB-8-1	Size 1	8/PK, 10 PK/BX
PKB-8-2	Size 2	8/PK, 10 PK/BX
PKB-8-3	Size 3	8/PK, 10 PK/BX
PKB-8-4	Size 4	8/PK, 10 PK/BX
PKB-8-5	Size 5	8/PK, 10 PK/BX
PKB-8-6	Size 6	8/PK, 10 PK/BX
PKB-8-A*	Assorted (Sizes 00-6)	8/PK, 10 PK/BX

***Assorted Pack:** A single, re-sealable packet contains 8 pieces – one brush for each size from 00 to 6



Interdental Brushes with Plastic Handles

- Perfect for cleaning in between teeth
- Suitable for daily at-home maintenance for healthy teeth and gums
- Plastic coated wire – safe around dental implants

Item#	Description	Packaging
PK-10-00	Size 00	10/PK, 10 PK/BX
PK-10-0	Size 0	10/PK, 10 PK/BX
PK-10-1	Size 1	10/PK, 10 PK/BX
PK-10-2	Size 2	10/PK, 10 PK/BX
PK-10-3	Size 3	10/PK, 10 PK/BX
PK-10-4	Size 4	10/PK, 10 PK/BX
PK-10-5	Size 5	10/PK, 10 PK/BX
PK-10-6	Size 6	10/PK, 10 PK/BX
PK-7-7	Size 7	7/PK, 4 PK/BX
PK-7-8	Size 8	7/PK, 4 PK/BX
PK-7-9	Size 9	7/PK, 7 PK/BX
PK-8-A*	Assorted (Sizes 00-6)	8/PK, 10 PK,BX
PK-4V-75**	Assorted (Sizes 00, 1, 3, 5)	




#PK-7-8
(7 pieces per pack,
all the same size)

***Assorted 8 Pack:** Each packet contains 8 pieces of the same size in re-sealable packet

****Assorted 4 Pack:** Each packet contains 4 brushes to cover most common used brush sizes, one brush for each size of 00,1, 3 & 5

New Products – Piksters Interdental Brushes and Protective Eyewear Interdental Brushes



Interdental Brushes with Plastic Handles

- Perfect for cleaning in between teeth
- Suitable for daily at-home maintenance for healthy teeth and gums
- Plastic coated wire – safe around dental implants

Item#	Description	Packaging
PK-40-00	Size 00	40/BX, 10 BX/CS
PK-40-1	Size 1	40/BX, 10 BX/CS
PK-40-2	Size 2	40/BX, 10 BX/CS
PK-40-3	Size 3	40/BX, 10 BX/CS
PK-40-4	Size 4	40/BX, 10 BX/CS
PK-40-5	Size 5	40/BX, 10 BX/CS
PK-40-6	Size 6	40/BX, 10 BX/CS
PK-40-7	Size 7	40/BX, 10 BX/CS



Professional Starter Kit

- Contains 116 brushes in a mix of sizes 00 to 9, size 1 Kink and SupaGrips. Used by the dentist to demonstrate to patients and determine the correct size.

Item#	Description	Packaging
PK-PSP-116	Starter Kit	116/BX



foamies™ Protective Eyewear

- Tinted and non-tinted disposable protective eyewear

Item#	Description	Packaging
UPK-F050C-S	Clear Lense - Small	50/BX, 10 BX/CS
UPK-F050C-M	Clear Lense - Medium	50/BX, 10 BX/CS
UPK-F050C-L	Clear Lense - Large	50/BX, 10 BX/CS
UPK-F050T-S	Tinted Lense - Small	50/BX, 10 BX/CS
UPK-F050T-M	Tinted Lense - Medium	50/BX, 10 BX/CS
UPK-F050T-L	Tinted Lense - Large	50/BX, 10 BX/CS



Bamboo Classic Toothbrushes

- 100% Biodegradable Bamboo Handle, each Toothbrush is packed in an individual package

Item#	Description	Packaging
PKB-CT	Bamboo Soft Bristle Toothbrush	10/BX, 4 BX/CS

Preventative- Prophy Angles & Fluoride Trays



Disposable Prophy Angles

- Precision gears for reduced vibration and fatigue
- Ribbed housing for increased maneuverability
- Smooth running and no cup dropping
- Not made with natural rubber latex
- Convenient 100-piece countertop dispenser box
- Available with firm or soft cups
- Non-sterile
- Single use

Packaging: 100 PC/DSPNSR, 5 DSPNSR/CTN, 500 PCS TOTAL, 4 CTN/CS

Item#	Description
UPC-9020-S	Prophy Angle (Green Cup) - Soft
UPC-9021-F	Prophy Angle (White Cup) - Firm



Fluoride Trays

- Hinged dual arch trays are color coded by size
- Locking flip-tab handle for easy filling, insertion and removal
- Strong, soft foam material for patient comfort and safety
- Available in sizes: small, medium, large and child/pedo
- Non-sterile
- Single use

Packaging: 100/BG, 10 BG/CS

Item#	Description
UFT-7401	Large - Blue
UFT-7402	Medium - White
UFT-7403	Small - Yellow
UFT-7404	Child/Pedo Size - Green



Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Find our full line of products at dukal.com