



Impregnated Gauze Dressings

Impregnated wound dressings are an important part of traditional wound care. These products maintain a moist wound environment to promote healing or enhance wound drainage. Dukal's line of impregnated wound dressings include oil emulsion dressings, petrolatum gauze, and petrolatum gauze with xeroform. These dressings require a secondary dressing, like our gauze sponges, non-woven dressings, or non-adherent pads. For wounds with depth, Dukal's packing strips are available in both plain and iodoform.

> For more information, please email customerservice@dukal.com

Petrolatum Dressings

- Provides soothing non-adherent protection with more convenience and less mess
- Promotes healing by maintaining a moist wound environment
- Clings and conforms to all body contours
- Reduces pain and trauma during dressing changes
- Decreases chances of skin maceration
- Available with Xeroform, a special deodorizing agent containing 3% bismuth tribomphenate



Item #	Description	Packaging
204-12	1/2" x 72"	12/BX
205-12	1" x 36"	12/BX
207-12	3" x 18"	12/BX
208-12	3" x 36"	12/BX
209-12	6" x 36"	12/BX
206	3" x 9"	50/BX
210	1" x 8"	50/BX
212	Xeroform, 1" x 8"	50/BX
213	Xeroform, 5" x 9"	50/BX
214	Xeroform, 2" x 2"	25/BX
215	Xeroform, 4" x 4"	25/BX
216	Xeroform, 4" x 3 yd	6/BX

Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com



Packing Strips

- Iodoform provides odor-resistance to the packing strip
- Tamper evident seals
- Not made with natural rubber latex

Item #	Description	Packaging
260	1/4" x 5 yd	12 BTL/CS
261	1/2" x 5 yd	12 BTL/CS
262	1" x 5 yd	12 BTL/CS
263	2" x 5 yd	12 BTL/CS
270	Iodoform, 1/4" x 5 yd	12 BTL/CS
271	Iodoform, 1/2" x 5 yd	12 BTL/CS
272	Iodoform, 1" x 5 yd	12 BTL/CS
273	Iodoform, 2" x 5 yd	12 BTL/CS



Oil Emulsion Dressings

- Made of a permeable knitted fabric that allows exudates to readily pass through to the second absorbent layer without adhering to granulation tissue
- Enhances wound drainage, supports a wound healing environment, and minimizes maceration
- Flexible and can be cut to size without tearing or unraveling

Item #	Description	Packaging
220	3" x 3"	50/BX
221	3" x 8"	24/BX
222	3" x 8"	36/BX
223	3" x 16"	36/BX
224	5" x 9"	12/BX



Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com