



Chemotherapy Gowns

Made of poly-coated SMS material to provide enhanced protection

- Hook and loop neck closure provides adjustability with a secure fit
- Waist ties and knit cuffs provide maximum comfort
- Our gown fabric was tested to ASTM F739-12 and passed for the below 16 chemotherapy drugs. No breakthrough was observed during test duration of 480 minutes based on average standardized breakthrough time:
 - Carmustine
 - Cisplatin
 - Cyclophosphamide
 - Cytarabine
 - Dacarbazine
 - Daunorubicin Hydrochloride
 - Doxorubicin Hydrochloride
 - Etoposide
 - Fluorouracil
 - Ifosfamide
 - Methotrexate
 - Mitomycin C
 - Mitoxantrone
 - Paclitaxel
 - Thiotepa
 - Vincristine Sulfate

Item #	Description	Packaging
G1005	Poly-Coated Chemotherapy Gown, Medium	10/BG, 3 BG/CS
G1006	Poly-Coated Chemotherapy Gown, Large	10/BG, 3 BG/CS
G1007	Poly-Coated Chemotherapy Gown, X-Large	10/BG, 3 BG/CS
G1008	Poly-Coated Chemotherapy Gown, XX-Large	10/BG, 3 BG/CS



Find our full line of products at [dukal.com](https://www.dukal.com)

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com