

Personal Protection



- **Quality**
- **Value**
- **Performance**





Gowns:

All gowns are made of high quality spunbonded polypropylene providing the best in basic coverage. Each gown features elastic cuffs for a more secure fit under gloves, and waist and neck ties for closure. Our Isolation Gowns are available in three different colors to allow easy detection between departments or functions. The Impervious Gowns feature a polyethylene coating making them fluid resistant providing additional protection. Available in yellow and blue. Our Patient Exam Gowns are open back and available in a deluxe version with waist tie. All DUKAL gowns are non-sterile and latex free.



301	Isolation Gown, Yellow	10/bg, 5 bg/cs
301WH	Isolation Gown, White	10/bg, 5 bg/cs
301BL	Isolation Gown, Blue	10/bg, 5 bg/cs
301SP	Isolation Gown, Yellow - Economy	10/bg, 5 bg/cs
303	Impervious Isolation Gown, Yellow	10/bg, 5 bg/cs
303BL	Impervious Isolation Gown, Blue	10/bg, 5 bg/cs
303W	Infection Control Gown Full Back, white	10/pk, 5 pk/cs
304	Patient Exam Gown, Dark Blue	10/bg, 5 bg/cs
305	Patient Exam Gown Deluxe, Dark Blue	10/bg, 5 bg/cs



Caps: Surgeon & Bouffant:

The classic Surgeon's Cap is made of a soft, reinforced material that is absorbent and comfortable to wear. The Nurse's Cap is made of a comfortable spunbonded polypropylene material which is cool and can be worn well for hours. All caps are latex free.



311	Surgeon's Cap, Blue	100/bx, 5 bx/cs
321B	Bouffant Cap, 21", Blue	100/bx, 5 bx/cs
324B	Bouffant Cap, 24", Blue	100/bx, 5 bx/cs
321W	Bouffant Cap 21", White	100/bx, 5 bx/cs
324W	Bouffant Cap, 24", White	100/bx, 5 bx/cs
321G	Bouffant Cap 21", Green	100/bx, 5 bx/cs
324G	Bouffant Cap, 24", Green	100/bx, 5 bx/cs



Disposable Scrub Wear & Coveralls:

DUKAL Scrubwear & Coveralls are made of a durable spunbonded polypropylene making them fluid resistant. All scrub pants have a drawstring tie waist, one front pocket and one back pocket and shirts have two pockets. They provide a great alternative to laundered scrubs. Our Coveralls are available in 3 sizes, have front zip closure, and elastic wrist and ankle cuffs. All DUKAL Scrubwear and Coveralls are latex free.

375S	Shirt, Small, Dark Blue	10/bg, 5 bg/cs
375M	Shirt, Medium, Dark Blue	10/bg, 5 bg/cs
375L	Shirt, Large, Dark Blue	10/bg, 5 bg/cs
375XL	Shirt, X-Large, Dark Blue	10/bg, 5 bg/cs
380S	Pants, Small, Dark Blue	10/bg, 5 bg/cs
380M	Pants, Medium, Dark Blue	10/bg, 5 bg/cs
380L	Pants, Large, Dark Blue	10/bg, 5 bg/cs
380XL	Pants, X-Large, Dark Blue	10/bg, 5 bg/cs
360	Shorts, Unisize, Dark Blue	10/bg, 5 bg/cs
382XL	Coveralls, White, XL	5/bg, 5 bg/cs
382XXL	Coveralls, White, XXL	5/bg, 5 bg/cs
382XXXL	Coveralls, White, XXXL	5/bg, 5 bg/cs
392XL	Coveralls, Blue, XL	5/bg, 5 bg/cs



Quality, Value & Performance is our PROMISE



Shoe Covers:

The Shoe Covers are made of a durable spunbonded polypropylene which makes them fluid repellent and rugged. They are latex free and a non-skid option is available.

350	Shoe Cover, Non-Skid, Blue	100/bx, 3 bx/cs
350E	Shoe Cover, Blue - Economy	100/bx, 3 bx/cs
352	Shoe Cover, Non-Skid, Extra Large, Size 14-16, Blue	100/bx, 2 bx/cs



Masks:

Procedure Masks:

DUKAL 3ply pleated Procedure Masks are made with lightweight breathable material and have a BFE of >97%. The Procedure Masks are available in earloop and tie on, with a flexible nosepiece for maximum comfort and fit. DUKAL masks are Latex Free.

Surgical Masks ASTM Moderate Protection*:

DUKAL 3 ply pleated Surgical Masks are made with lightweight breathable material, and have a BFE of >99%, and meet or exceed ASTM1 Standards. The Surgical Masks are available in earloop and tie on, with a flexible nosepiece for maximum comfort and fit. DUKAL masks are Latex Free.


N95 Masks ASTM 1 High Protection*:

The DUKAL Respirator/ Surgical Masks meet NIOSH and CDC guidelines. As a respirator, they are designed to reduce wearer exposure to airborne diseases and meet the CDC guidelines for tuberculosis exposure control with a filter efficiency of 95% or greater against particulate aerosols free of oil (Type N95 Respirator).

As surgical masks they are designed to provide a physical barrier to fluids and particulate materials with a High Fluid Penetration Resistance (ASTM F1862-05 / 160mmHg), a Particle Filtration Efficiency (PFE) of >99% and a Bacterial Filtration Efficiency (BFE) of 99.9%.

Our highly effective particulate respirator/surgical masks are available in duckbill (folded), cone style, and cone with a one way valve.

1530	Procedure Mask w/Ear Loop, Blue	50/bx, 6 bx/cs
1531	Procedure Mask w/Ear Loop, Yellow	50/bx, 6 bx/cs
1540	Surgical Mask w/Tie, Blue	50/bx, 6 bx/cs
1541	Surgical Mask w/Ear Loop, Blue	50/bx, 6bx/cs
1542	Surgical Mask w/Ear Loop, Yellow	50/bx, 6bx/cs
1550	Anti-Fog Surgical Mask w/Tie, Blue	50/bx, 6 bx/cs
1560	Surgical Mask w/Fluid Shield, w/Ear Loop, Blue	25/bx, 8 bx/cs
1570	Folded, White, N95	20/bx, 20 bx/cs
1571	Cone, White, N95	20/bx, 20 bx/cs
1572	Cone, One-Way Valve, White	5/bx, 20 bx/cs



Product Categories

Gauze Products

Gauze Sponges
Basic Gauze Sponges
Super Sponges
Gauze Pads
Tray Gauze
USP Type VII
X-Ray Detectable Gauze USP Type VII
Cotton Filled Gauze

Non-Woven Dressing

Clinisorb Non-Woven Sponges
New Sponges
Premium Sponges
I.V. & T-Drain Sponges

Specialty Dressing

Top Sponges,
Emergency Pressure Dressing
Combine ABD Pads
Non-Adherent Pads
Eye Pads
Trauma Dressing
Post-Op Sponges
Pressure Dots

Advanced Wound Care

Wound Closure Strips
Transparent Semi Permeable Dressing
Suture Removal Kits
Tracheostomy Care Kits
Take Home Wound Care Kits
Caliber™ Non Woven Island Dressings
Caliber™ Oil Emulsion Dressings

Bandages

Rolled Gauze
Elastic Bandage Roll
Conforming Stretch Gauze
Basic Conforming Stretch Gauze
Caliber™ Adhesive Bandages
Cohesive Bandages
Bandage Roll
Triangular Bandages

Cotton & Wood Products

Cotton Stockinette
Cotton Tipped Applicators
Cotton Balls
1 lb Cotton Roll
Tongue Depressors
Wood Applicators

Tapes

Cloth Tape
Paper Tape
Transparent Tape
Waterproof Tape
Single Use Tapes

Specialty Surgical Dressing

Cherry Sponge
Kittner Dissectors
Laprotomy Sponges
OR Towels
Peanut Sponges
Tonsil Sponges
Flat Stick Sponges
Peanut Sponges
Cherry Sponges
Caliber™ Surgical Drapes
Burn Dressings

Personal Protection

Caps
Gowns
Masks
Shoe Covers
Disposable Scrub Wear & Coveralls

Disposable Linens

Disposable Linens
Emergency Blankets
Survival Wrap
Burn Sheet

Hygiene

Alcohol Prep Pads-Medium
Dry Wash Cloths
Fluff Underpad
Wet Wipes
Disposable Vaginal Speculums
BZK Swabs
BZK Towelettes
Cleansing Towelettes
Sting Relief Pads

Health & Beauty

Atmosphere™ Breathing Extender
DUKAL Reflections™ Chem Kit
DUKAL Reflections™ Wood Applicators
DUKAL Reflections™ Glycolic Applicators
DUKAL Reflections™ Esthetic Wipes
Body Toaster™ Spa Wrap
DUKAL Reflections™ Disposable Undergarments

