



PASSION. PARTNERSHIP. POSSIBILITIES.



# Caliber™ Adhesive Bandages

## Description

Durable plastic, sheer, or flexible fabric bandages to meet all market needs

## Features & Benefits

- Absorbent, non-adherent pad protects and cushions the wound without sticking
- Sterile and individually wrapped in easy to peel pouches
- Not made with natural rubber latex

| Item#          | Description   | Packaging        | HCPCS Code |
|----------------|---|------------------|------------|
| <b>Sterile</b> |   |                  |            |
| 7604           | Sheer, 3/4" x 3"  | 100/bx, 24 bx/cs | A6413      |
| 7608           | Sheer, 1" x 3"  | 100/bx, 24 bx/cs | A6413      |
| 7612           | Sheer, 2" x 4", X-Large                                 | 100/bx, 24 bx/cs | A6413      |
| 7615           | Sheer, Spot, 7/8" Round                                 | 100/bx, 24 bx/cs | A6413      |
| 7605           | Flexible Fabric, 3/4" x 3"                              | 100/bx, 24 bx/cs | A6413      |
| 7609           | Flexible Fabric, 1" x 3"                                | 100/bx, 24 bx/cs | A6413      |
| 7613           | Flexible Fabric, 2" x 4", X-Large                       | 100/bx, 24 bx/cs | A6413      |
| 7618           | Flexible Fabric, Fingertip                              | 100/bx, 24 bx/cs | A6413      |
| 7602           | Flexible Fabric, Knuckle                                | 100/bx, 24 bx/cs | A6413      |
| 7616           | Plastic, 3/4" x 3"                                      | 100/bx, 24 bx/cs | A6413      |
| 7617           | Plastic, 1" x 3"  | 100/bx, 24 bx/cs | A6413      |
| 7611           | Plastic, 2" x 4", X-Large                               | 100/bx, 24 bx/cs | A6413      |
| 7614           | Plastic, Spot, 7/8" Round                               | 100/bx, 24 bx/cs | A6413      |
| 7600           | Butterfly Wound Closure Strips, Medium, 3/8" x 1-13/16" | 100/bx, 24 bx/cs | A4450      |

### Connect with us!

Find our full line of products online at:  
[www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
 Toll Free: 1-800-243-0741  
 Fax: 1-631-656-3810  
 Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
 To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
 2 Fleetwood Court,  
 Ronkonkoma, NY 11779

