



# Surgical Drapes



**DUKAL's 3 ply Caliber™ Surgical Drapes are excellent for use during minor surgical procedures. Two layers of absorbent tissue with a polyethylene lining in between help eliminate risk of bacterial migration and moisture strike through. Available in both plain for general use, and fenestrated to allow access to operative site while protecting the surrounding area. Sterile, latex free.**

ITEM #	DESCRIPTION	PACKAGING
20-001	Sterile, Fenestrated Surgical Drape, 18" x 26", 3ply	1/pk, 50/bx, 6bx/cs
20-002	Sterile, Surgical Drape, 18" x 26", 3ply	1/pk, 50/bx, 6bx/cs

Quality



Value



Performance

