



PASSION. PARTNERSHIP. POSSIBILITIES.



Retention Tape

Description

An adhesive non-woven fabric tape that is highly conformable

Features & Benefits

- Soft & flexible making it ideal for difficult to tape areas
- Printed s-curve release liner
- Pre-cut to eliminate waste
- Provides light compression without constriction

Suggested Uses

- To secure primary dressings
- To secure catheters and drainage tubes

Item	Description	Packaging
Non-Sterile - Latex Free		
8220	2" x 11 yds	1 rl/bx, 10 bx/cs
8221	4" x 11 yds	1 rl/bx, 6 bx/cs
8222	6" x 11 yds	1 rl/bx, 4 bx/cs

Connect with us!

Find our full line of products online at:
www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779

