



PASSION. PARTNERSHIP. POSSIBILITIES.



Basic Care Bandage Rolls

Description

Low ply, lightweight alternative to our standard Fluff Bandage Rolls



Features & Benefits

- Made of washed, fluff-dried 100% cotton gauze
- Open, crinkle weave provides quick wicking, maximum absorbency, and greater aeration
- Finished edges reduces linting and loose threads
- Adds bulk and loft that will conform as needed, secure dressings in place, and cushion the wound



Item	Description	Packaging	HCPCS Code
Sterile - Latex Free			
* 8642	2" x 4.5 yds, 3-ply	1 r/bg, 96 r/cs	A6445
* 8643	3" x 4.5 yds, 3-ply	1 r/bg, 96 r/cs	A6446
8645	4.5" x 4.1 yds, 3-ply	1 r/bg, 100 r/cs	A6446
* 8644	6" x 4 yds, 3-ply	1 r/bg, 48 r/cs	A6447

* New sizes

Connect with us!

Find our full line of products online at:
www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

