



PASSION. PARTNERSHIP. POSSIBILITIES.



Medicine Cups

Features & Benefits

- A convenient way to distribute medicine in pill or liquid format
- Calibrated with 5 measurements – ounces, drams, cubic centimeters, milliliters, and tablespoons
- Easy to read gradations
- Made of unbreakable translucent polypropylene



Item	Description	Packaging
9007	1 oz., Translucent	100/bg, 50 bgs/cs

Connect with us!

Find our full line of products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

