



PASSION. PARTNERSHIP. POSSIBILITIES.



albahhealth Xeroform Petrolatum Gauze

Description

Provides soothing non-adherent protection with more convenience, less mess, and freshens as it protects

Features & Benefits

- Non-adherent petrolatum
- Clings and conforms to all body contours
- Non-toxic, non-irritating, guaranteed sterile
- Contains 3% bismuth tribomophenate, a special deodorizing agent
- Not made with natural rubber latex

Common Applications:

- An initial layer in surgical dressings
- Skin graft sites
- Treatment of first degree burns
- Abrasions
- Wound packing

| Item# | Description | Packaging |
|----------------|-------------|-----------|
| Sterile | | |
| 212 | 1" x 8" | 50/bx |
| 213 | 5" x 9" | 50/bx |
| 214 | 2" x 2" | 25/bx |
| 215 | 4" x 4" | 25/bx |
| 216 | 4" x 3 yds | 6/bx |

Connect with us!

Find our full line of products online at: www.dukal.com

Telephone: 1-631-656-3800
Toll Free: 1-800-243-0741
Fax: 1-631-656-3810
Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779

