

Tech-Med[®] Single Post Instrument Stand

Features & Benefits

- Includes highest grade stainless steel tray (18/8,304)
- U-shaped base minimizes interference with other equipment
- Secure height adjustments 28"-48"
- Tray Size: 12 1/2" x 19"
- Base Width: 19", depth: 22"
- Weight Capacity: 11 lbs
- Replacement tray available



Item#	Description	Packaging
Non-Sterile		
4352	Single Post	1/ea
4264	Replacement Tray, 12-1/2" x 19"	1/ea

Tech-Med[®] Mobile Instrument Stand

Features & Benefits

- Includes highest grade stainless steel tray (18/8, 304)
- Mobile base features locking casters
- Secure height adjustments 28" – 48"
- Tray size: 12-1/2" x 19"
- Base width: 20"
- Weight Capacity Standard: 11 lbs
- Replacement tray available



Item#	Description	Packaging
Non-Sterile		
4365	Mobile Base	1/ea
4264	Replacement Tray, 12-1/2" x 19"	1/ea

Connect with us!

Find our full line of products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779



Tech-Med[®] Double Post Mayo Stand with Large Tray

Features & Benefits

- Includes highest grade stainless steel tray (18/8, 304) that is 25% larger than standard trays and ideal for urgent care facilities
- Base minimizes interference with other equipment
- Secure height adjustments 37" – 53"
- Tray size: 16-3/4" x 21-1/2"
- Base width: 21", depth: 25"
- Weight Capacity Standard: 50 lbs
- Replacement tray available



Item#	Description	Packaging
Non-Sterile		
4366	Double Post	1/ea
4266	Replacement Tray, 16-3/4" x 21-1/2"	1/ea

Tech-Med[®] Double Post Mayo Stand

Features & Benefits

- Includes highest grade stainless steel tray (18/8, 304)
- Base minimizes interference with other equipment
- Secure height adjustments 34" – 53"
- Tray size: 12-1/2" x 19"
- Base width: 20", depth: 19-1/2"
- Weight Capacity Standard: 15 lbs
- Replacement tray available



Item#	Description	Packaging
Non-Sterile		
4368	Double Post	1/ea
4264	Replacement Tray, 12-1/2" x 19"	1/ea

Connect with us!

Find our full line of products online at:
www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779

