


# Advazorb Fixation<sup>®</sup>

## Tracheostomy Dressing



## Gentle on skin

Distributed by



PASSION. PARTNERSHIP. POSSIBILITIES.

[www.dukal.com](http://www.dukal.com)





[www.advancismedical.com](http://www.advancismedical.com)

# Advazorb Fixation<sup>®</sup> Tracheostomy Dressing

Advazorb Fixation is a post decannulation tracheostomy dressing fitted with a central location button.

Advazorb Fixation is a hydrophilic foam dressing with a soft silicone contact layer and surrounding adhesive. It is a low profile, non-bulky, skin friendly foam dressing.

- The silicone adhesive combines security with pain-free dressing changes. The waterproof backing enables the patient to shower without the need to replace or cover the dressing.
- The central location button prevents unwanted escape of air when gentle pressure is applied during talking, laughing, coughing and sneezing.



## Advazorb Fixation<sup>®</sup> Tracheostomy Dressing

The silicone wound contact layer has large pores which enable the passage moisture into the absorbent foam while protecting the healing stoma site. This combination ensures the dressing stays comfortably in place while minimizing peristomal irritation.

- Supplied sterile, for single patient use
- Change every 5-7 days, or as needed
- Allows atraumatic dressing change
- Reduces medical adhesive-related skin injury (MARSI)
- Ease of use creates consistency upon application
- Provides a barrier for bacteria
- Latex free

Item #	Description	Packaging
CR4341	4.5 in x 4 in	10/bx

Manufactured in the UK, the Advancis Medical<sup>®</sup> range of Advazorb<sup>®</sup> Hydrophilic foam dressings, offer the combination of a high quality performance product at a competitive price.



DUKAL Corporation  
2 Fleetwood Court  
Ronkonkoma, NY 11779  
Tel: 1-800-243-0741 | Email: [advancisinfo@dukal.com](mailto:advancisinfo@dukal.com)

**Advancis**  
Medical  
[www.advancismedical.com](http://www.advancismedical.com)