



PASSION. PARTNERSHIP. POSSIBILITIES.



# Silver Nitrate Applicators

## Description

Used for cauterization of skin or mucus membrane and for the removal of granulation tissue, warts and verrucae

## Features & Benefits

- 75% silver nitrate, 25% Potassium Nitrate
- Plastic caustic applicator may be bent or shaped to enable easy access to the ear, nose or throat cavities
- The flexible 6" plastic caustic applicator is molded at high temperatures which make the applicator hygienic and splinter free
- Federal law prohibits dispensing without prescription
- Five year shelf life

## Common Applications:

- Nose Bleeds
- Vasomotor Rhinitis
- Post Bioptic Cervical Nicks

| Item# | Description | Packaging             |
|-------|-------------|-----------------------|
| 7482  | 6"          | 100/vial, 10 vials/bx |

### Connect with us!

Find our full line of products online at: [www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
Toll Free: 1-800-243-0741  
Fax: 1-631-656-3810  
Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
2 Fleetwood Court,  
Ronkonkoma, NY 11779



\*Bray Healthcare is a registered trademark of Bray Group Ltd.